ALT - TGP

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
ENGLISH
INTENDED USE
Kit intended to determine ALT/GPT activity in serum and plasma (EDTA and heparin). Diagnostic use only.
STORAGE AND HANDLING
• Store at 2 to 8°C
• Protect from light.
• Stable until the kit expiration date that is printed on the label.
• Do not use reagents whose date has expired.

WORKING PRINCIPLE

The Alanine Transaminase (ALT or GPT) catalyses the transference of amino group from alanine to alpha-ketoglutarate. This transference produces pyruvate and glutamate. The activity is then determined by coupling this reaction to Lactate Dehydrogenase (LDH) enzyme that transform the pyruvate product in lactate consuming NADH which is measured at 340 nm

SAMPLE PREPARATION AND STABILITY

Sample Type: Serum and plasma (EDTA and heparin)
Collection, handling and preparation: Make the collection of the sample in accordance with Good Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: The ALT-TGP in serum and plasma is stable for 7 days if stored at a temperature of 4 to 8°C or -20°C

PRODUCT DESCRIPTION

R 1	Tris buffer ≥ 100 mM L-Alanine $> 0,5$ M Lactate Dehydrogenase ≥ 1500 U/L; triton X100 $\geq 0,02\%$ alpha-ketoglutarate ≥ 10 mM activators, stabilizers, preservative.	
R 2	Carbonat buffer ≥ 1 mM NADH $\geq 0,5$ mM preservative.	

The determination of ALT-TGP is traceable to reference material ERM AD454/IFCC

STABILITY IN USE

- The stability of the product (R1 and R2) in use is 18 months as long as followed by the recommended storage conditions (2 to 8°C).
- The stability of working reagent is 14 days, as long as followed by the conditions of preparation and recommended storage (2 to 8°C).
- Reagents must remain outside the specified temperature only the time required for testing

TECHNICAL PREPARATION

A) REAGENT PREPARATION

Work Reagent (WR)

Mix in proportion of 4 parts of R1 + 1 part of R2.

Homogenize gently. The Work Reagent is stable 14 days at 2 to 8°C

B) OPERATING RANGE

The product operating range is from 5,15 U/L to 350 U/L. For higher values, dilute the sample with NaCl 150 mM (0,9%), repeat the assay and multiply the result by the dilution factor.

QUALITY CONTROL

Use of controls should be a routine practice in the lab. It is suggested to use a control within the reference range or decision level and another control with value in another range of clinical significance. For the internal quality control of the lab it's indicated the use of the calibrator serum and control sera below

Calibrador Serum - Autocal H	REF	13.002.00
Normal Control Serum - Quantinorm		13.003.00
Pathological Control Serum - Quantilat		13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Heat the reagent for 3 minutes at 37°C
2. Pipette in the assay tube:

Work Reagent	Volume
Sample	1,0 mL
	100 µL

3. Homogenize and insert immediately in the thermostatic cuvette at 37°C
4. After 1 minute note the initial absorbance (A_0) and read again each minute, during 3 minutes (A_1, A_2 and A_3 , respectively).

B) CALCULATION

Using the measured absorbances, calculate the mean variation of absorbance per minute ($\Delta A/\text{min}$):

$$\Delta A/\text{min} = \frac{(A_0 - A_1) + (A_1 - A_2) + (A_2 - A_3)}{3}$$

The activity of the GPT is then calculated by multiplying the $\Delta A/\text{min}$ by the factor below

$$GPT (U/L) = \Delta A/\text{min} \times 1746$$

Example:

$A_0 = 1,297$; $A_1 = 1,278$; $A_2 = 1,263$; $A_3 = 1,245$

$$\Delta A/\text{min} = \frac{(1,297 - 1,278) + (1,278 - 1,263) + (1,263 - 1,245)}{3}$$

$$GPT (U/L) = 0,017 \times 1746 = 29,7 U/L$$

C) INTERPRETATION

The Glutamate Pyruvate Transaminase (GPT) or Alanine Transaminase (ALT) is an intracellular enzyme present in high quantities in liver, kidney and in low quantities in skeletal and cardiac muscles. As a liver function marker, the GPT is more sensitive to the

detection of hepatocyte injury than to gall bladder obstruction and is considered an excellent cellular marker. Elevated levels of GPT are found in infections and intoxication of the liver, pancreatic diseases, cirrhosis, obstructive jaundice and metastasis. In patients with myocardial infarction the GPT is normal or slightly elevated. However, cardiac insufficiency with consequent liver necrosis may lead to elevated values of GPT.

INTERFERING AND LIMITATIONS

Anti-coagulants: Citrate, sodium fluoride and sodium oxalate may interfere.

Hemolysis, Jaundice and Lipemic Sera: Hemoglobin > 200 mg/dL, Bilirubin > 18 mg/dL, Triglyceride > 750 mg/dL.

PERFORMANCE CHARACTERISTICS

Sensitivity: Detection limit: 2.69 U/L / quantification limit: 5.15 U/L

Analytical Specificity: The product determines ALT-TGP specifically in the presence of other interfering substances in the sample up to the concentrations stated above.

Accuracy: The method was compared with a similar method for the determination of 40 samples in duplicate. The regression equation obtained was $y = 1.006x - 0.53$ and the correlation coefficient $r = 0.9993$. Using this equation on the total systematic error estimates -0.72% to a level of 40 U/L and 0.07% to a level of 100 U/L.

Precision: Determined using samples in 03 levels of decision, with two runs in duplicate per day during 20 days, being obtained:

Samples (U/L)	Repetitions	Within-Run Precision		Total Precision	
		SD	%CV	SD	%CV
22,80	80	0,48	2,1	0,59	2,6
47,20	80	0,73	1,5	1,18	2,5
123,12	80	1,14	0,9	1,62	1,3

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK WARNINGS AND PRECAUTIONS

- Use of Protective Equipment in accordance with Good Laboratory Practice
- Do not mix reagents from different lots.
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any significant disagreement with the one specified in the product MSDs.
- Avoid leaving reagents outside the specified storage conditions, when they are not in use.
- Use glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The level of the water bath must be greater than the test tubes containing the reaction.
- The level of water in the thermostat bath must be over the level of mix reaction with in the assay tubes.
- The laboratory should establish chemical, microbiological and particulate requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

REFERENCE RANGES

NORMAL VALUES	37 °C
ME N	< 42 U/L
WOMEN	< 32 U/L

These values are intended for orientation only. It's recommended that each lab establishes its own reference range.

Conversion for International System of Units (SI): GPT (U/L) X 0,017 = GPT (µkat/L)

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer with a thermostated cuvette able to read at 340 nm
- Thermostatic water bath at 37°C
- Glass pipettes and/or automatic
- Clock or chronometer.
- Assay tubes

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.bioteccalind.com or calling for 55-35-3214-4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use Biotécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiry date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@bioteccaltda.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.bioteccaltda.com.br

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- SOCIETÀ CHIMICA ITALIANA - WORKING GROUP ON ENZYMES. International Federation of Clinical Chemistry/IFCC methods for measurement of catalytic concentration of enzymes. **Clin Chim Acta** v. 281, p.1-2, p.55-539, 1999.
- BURTONS, C A; ASHWOOD, E R. Tietz textbook of clinical chemistry, 4 ed. Philadelphia: W B Saunders, 1994, 836p.
- YOUNG, D.S. **Effects of drugs on clinical laboratory tests - v. 2**, 5 ed. Washington DC: AACCPress, 2000.
- WESTGARD, J. O et al. A multi-rule shewhart chart quality control in clinical chemistry. **Clin Chem** v.27 p.493-501, 1981.

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