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
ENGLISH
INTENDED USE
Kit intended to determination of uric acid in serum and urine. Diagnostic use only.

- STORAGE AND HANDLING**
- Store at 2 to 8 ° C.
 - Protect from light.
 - Stable until the kit expiration date that is printed on the label.
 - Do not use reagents whose date has expired.

WORKING PRINCIPLE

The uric acid from sample in presence of uricase and dioxygen produces allantoin and hydrogen peroxide. The hydrogen peroxide, in presence of DHBS and 4-aminoantipyrine, reacts with Peroxidase producing a pink compound (Quinonimine) that absorbs the light at 505 nm. The color intensity is proportional to the concentration of uric acid in the sample.

SAMPLE – PREPARATION AND STABILITY
Sample Type: Serum and urine.
Collection, handling and preparation: Make the collection of the sample in accordance with Good Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation:
Serum: The uric acid is stable in serum for 5 days from 4 to 8 °C and 6 months from -20 °C. The uric acid in urine 24 hours is stable for 4 days from 20 to 25 °C.
Urine: Work with samples collected in the period of 24 hours. Do not add preservatives. Homogenize the urine, measure the volume and separate 20 mL. Adjust the pH between 7 and 9 with 5% NaOH solution and heat the sample for 10 minutes at 56 °C. Prepare a 1:10 dilution of the urine with purified water and perform the assay. Multiply the result obtained by 10. This dilution should result in a measurement within the operating range. For higher values, perform new dilution by changing the ratio.

PRODUCT DESCRIPTION

R 1	Phosphate buffer ≥ 50 mmol/L; Uricase ≥ 100 U/L; Peroxidase ≥ 1000 U/L; 4-aminoantipyrine ≥ 50 µmol/L; DHBS ≥ 0,5 mmol/L; Sodium Azide 0.05% w/v; activators; stabilizer and preservative.	X
STD	Uric acid concentration of 10 mg/dL. Isopropilic alcohol ≥ 10% v/v; stabilizer and preservative.	X

The determination of Uric Acid is traceable to reference material NIST 913a.

- STABILITY IN USE**
- The stability of the product (R1 e STD) in use is 14 months as long as followed by the recommended storage conditions (2 to 8 °C).
 - Reagents must remain outside the specified temperature only the time required for testing

TECHNICAL PROCEDURE
A) REAGENT PREPARATION
Reagent 1 (R1)
Ready-to-use Reagent
B) OPERATING RANGE
Serum: The product operating range is from 0,23 mg/dL to 20,00 mg/dL. For higher values, dilute the sample with NaCl 150 mM (0,9%), repeat the assay and multiply the result by the dilution factor.
Urine: The operating range of the product is 9,19 mg / dL to 210,10 mg / dL.

QUALITY CONTROL
Use of controls should be a routine practice in the lab. It is suggested to use a control within the reference range or decision level and another control with value in another range of clinical significance. For the internal quality control of the lab it's indicated the use of the calibrator serum and control sera below:
Calibrator Serum - Autocal H
Normal Control Serum – Quantinorm
Pathological Controle Serum – QuantiAlt

	REF	13.002.00
		13.003.00
		13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION
A) TEST PROCEDURE
1. Pipette in the assay tubes:

	Blank	Standard	Sample
STD	-	20 µL	-
Sample	-	-	20 µL
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate 10 minutes at 37 °C.
3. Measure the absorbance of Standard (Ast) and Sample (Asa) against the blank at 505 nm (490-510 nm). The final reaction is stable for 15 minutes.

B) CALCULATIONS
Calculations of serum
Uric Acid (mg/dL) = $\frac{\text{Sample Absorbance}}{\text{STD Absorbance}} \times \text{STD Concentration (mg/dL)}$

With Calibration Factor (CF):
CF = $\frac{\text{STD Concentration (mg/dL)}}{\text{STD Absorbance}}$
Uric Acid (mg/dL) = Sample Absorbance x CF

Calculations of Urine:
Urine (mg/24 hours) = $\frac{* \text{mg/dL} \times \text{urinary volume (in mL)}}{100}$

* Determine the concentration of uric acid in the urine sample in mg/dL using the same calculation procedure for serum. Multiply the result by 10.

C) INTERPRETATION
Uric acid is the end product of purine metabolism in human organism. Its dose is used predominantly in the gout, clinical condition which is established as a result of a primary or hyperuricemia caused by other conditions or treatments which cause secondary hyperuricemia. Furthermore, the determination of uric acid is used in the diagnosis and monitoring of numerous renal and metabolic disorders, including renal failure, leukemia, psoriasis, malnutrition, in patients using cytotoxic drugs, and the pregnancy-induced hypertension (preeclamptic toxemia).

INTERFERING AND LIMITATIONS

Anticoagulantes: Citrate, EDTA, fluoride and oxalate may interfere in the reaction.

Hemolyzed, Jaundice and Lipemic Sera:

Serum: Hemoglobin > 200 mg/dL, Bilirubin > 6,0 mg/dL, Triglyceride > 900 mg/dL.

Urine: Hemoglobin > 200 mg/dL / Bilirubin > 15 mg/dL.

Drugs: see recommended reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity:

Serum: Detection limit: 0,18 mg/dL / quantification limit: 0,23 mg/dL.

Urine: Detection limit: 0.146 mg / dL / Limit of quantification: 9,19 mg/dL

Analytical Specificity: The product determines uric acid specifically in the presence of other interfering substances in the sample, up to the concentrations stated above

Accuracy:

Serum: The method was compared with a similar method for the determination of 40 samples in duplicate. The regression equation obtained was $y = 0,991x + 0,007$ and the correlation coefficient $r = 0,9901$. Using this equation the total systematic error estimated is $-0,81\%$ to a level of 8,00mg/dL and $-0,85\%$ to a level of 15,00mg/dL.

Urine: The method was compared to a similar method by the determination of 27 samples in duplicate. The regression equation $y = 0,987x - 1,762$ and correlation coefficient $r = 0,9944$ was obtained. Using this equation, the total systematic error estimated was $-1,89\%$ for a level of 300,0 mg/24hs and $-1,57\%$ for a level of 650,0 mg/24h.

Precision:

Serum: Determined using samples in 03 levels of decision, with two runs in duplicate per day during 20 days, as follows:

Samples (U/L)	Repetitions	Within-Run Precision		Total Precision	
		SD	%CV	SD	%CV
3,891	80	0,026	0,7	0,057	1,5
6,391	80	0,054	0,8	0,090	1,4
12,506	80	0,070	0,6	0,185	1,5

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

Urine: It was determined using samples at 02 decision levels, two runs in duplicate per day for 20 days, as follows:

Samples (mg/dL)	Repetitions	Within-Run Precision		Total Precision	
		SD (mg/dL)	%CV	SD (mg/dL)	%CV
13,264	80	0,265	2,0	0,332	2,2
19,406	80	0,498	2,6	0,515	2,7

% CV: Coefficient of variation expressed as a percentage; SD: Standard deviation.

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with Good Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any sign in disagreement with the one specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions, when they are not in use.
- Use glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The level of the water bath must be greater than the test tubes containing the reaction.
- The level of water in the thermostatic bath must be over the level of mix reaction within the assay tubes.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

REFERENCE RANGES

	MEN		WOMEN	
	mg/dL	mmol/L	mg/dL	mmol/L
Serum	3,5 - 7,2	208 - 428	2,6 - 6,0	155 - 357
Urine	250 - 750 mg/24 hours			

These values are intended for orientation only. It's recommended that each lab establishes its own reference ranges.

Conversion for International System of Units (SI): mmol/L

Uric Acid (mg/dL) X 59,5 = mmol/L

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer with a thermostated cuvette able to read at 505 nm (490 - 510 nm).
- Thermostatic water bath at 37°C.
- Glass Pipettes and/or automatic.
- Clock or chronometer.
- Assay tubes

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotecnica.ind.br or calling for 55-35-3214-4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on

handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotecnicaltda.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.biotecnicaltda.ind.br

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