## BIURRUBI NA DRETA E TOTAL

Responsible Técnico: Dr. Gílson Sérgio Rozo
CRF MG - 5310
MS 80027310205

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ENGLISH
INTENDED USE
Kit for determination of Direct and Total Bilirubin in serum In vitro diagnostic use.

STORAGE AND HANDLING
- Store at 2 to 8 °C
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE
The determination is based on the azo reaction, which colorimetrically quantifies the reaction of azobilirrubina when bilirubin reacts under certain conditions with azotized sulfanilic acid. The intensity of the formed color is proportional to the concentration of bilirubin in the sample and in the presence of methyl sulfoxide (accelerator), the indirect form (unconjugated) and the direct form (nonconjugated, conjugated and protein bound) are reacted to determine Total Bilirubin. The direct form is determined without the accelerator.

SAMPLES - TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum
Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material. Immediately after collection, protect the sample from the action of light.
Preservation: Serum Bilirubin is stable for 7 days at 4 to 8 °C and 6 months at -20 °C.

Table with 2 columns: Symbol (RBD, RBT, RNI) and Description (Sulfamic acid > 10mmol/L, hydrochloric acid > 0.5% v/v, etc.)

STABILITY IN USE
- Once opened, the product (RBD, RBT and RN) in use is stable to the expiration date printed on the label, since recommended storage conditions (2 to 8 °C) are followed.
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

TREATMENT AND HANDLING OF THE PRODUCT

A) REAGENT PREPARATION
RBD, RBT e RN: Ready-to-use reagents.

B) OPERATIONAL INTERVAL
The operational range of Direct Bilirubin is from 0.17 mg / dL to 15.0 mg / dL. The operational range of Total Bilirubin is from 0.18 mg / dL to 15.0 mg / dL. For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL
The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For Calibration and Internal Laboratory Quality Control it is recommended to use the calibrator serum and the control sera below.
Serum Calibrator - Autocal H 13.002.00
Norma Control Serum Quantinorm 13.003.00
Serum Pathologic Control - Quanti d t 13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

Table with 4 columns: Symbol (RBD, RN, Sample/Calibrador) and Description (B ank, Blank Sample, Sample). Rows for 1,0 mL, 30 µL, 50 µL.

2. Homogenize well and keep the tubes for 5 minutes at room temperature.
3. Measure the absorbance of the Calibrador, Calibrador Blank, Sample and Sample Blank again at 550 nm (530-570). Color is stable for 5 minutes.

Table with 4 columns: Symbol (RBT, RN, Sample/Calibrador) and Description (B ank, Blank Sample, Sample). Rows for 1,0 mL, 30 µL, 50 µL.

2. Homogenize well and keep the tubes for 5 minutes at room temperature.
3. Measure the absorbance of the Calibrador, Calibrador White, Sample and Sample Blank again at 550 nm (530-570). Color is stable for 30 minutes.

C) CALIBRATION
Use Calibrator Autocal H Ref. 13.002.00. The concentration of Bilirubin in the calibrator is traceable to the Standard Reference Material (SRM) 916a of the National Institute of Standards and Technology (NIST).

D) CALCULATIONS
Bilirrubina (ng / dL) = (Sample Abs - Sample White Abs) x Calibrador Value (ng / dL) / (Calibrador Abs - White Calibrador Abs)

Exemplo
Calibrador Concentration = 4,62 mg / dL
Sample Absorbance = 0,194
Sample Blank Absorbance = 0,040
Calibrador Absorbance = 0,234
Absorbance Blank Calibrador = 0,024
Bilirrubina (ng / dL) = 0,194 - 0,040 x 4,62 = 3,39 ng / dL
0,234 - 0,024

With Calibrators Factor
Calibrators Factor = Calibrador Concentration (mg / dL) / (Calibrador Abs - White Calibrador Abs)

Bilirrubina (ng / dL) = (Sample Abs - Abs Sample White) x Calibrators Factor
Exemplo
Calibrators Factor = 4,62 = 22
0,234 - 0,024
Bilirrubina (ng / dL) = (0,194 - 0,040) x 22 = 3,39 ng / dL.

Q) INTERPRETATION
Bilirubin is a pigment derived from hemoglobin degradation. Bilirubin is transported to the liver, poorly associated with albumin, in its unconjugated form. Carried into the hepatocyte it is rapidly conjugated to produce the bilirubin glucuronides, which are then excreted into the bile. Elevations of conjugated (direct) bilirubin are usually specific markers of hepatic dysfunction.

acute hepatitis, or cholestasis. Increased unconjugated (indirect) bilirubin is usually caused by hemolytic processes or by a reduction in liver metabolism.

**INTERFERING AND LIMITATIONS**

**Hemolysis, Jaundice and Lipemia:** Hemoglobin > 150 ng / dL, hemolysis interferes in the dosages of direct bilirubin and total bilirubin.

**Medications:** consult recommended literature reference (Young, 2000).

**PERFORMANCE CHARACTERISTICS**

**Sensitivity**

	Direct Bilirubin	Total Bilirubin
Detect on Limit	0.14 mg/dL	0.13 mg/dL
Qualification Limit	0.17 mg/dL	0.18 mg/dL

**Analytical Specificity:** The product specifically determines Direct and Total Bilirubin in the presence of other interfering substances in the sample up to the concentrations reported above.

**Accuracy:** The method was compared with a similar method by the determination of 40 samples in duplicate. The results are as follows:

**Direct Bilirubin:** Regression equation  $y = 1.017x - 0.008$  and correlation coefficient  $r = 0.9909$ . Using this equation the estimated total systematic error of 0.9% for a level of 1.0 mg / dL and 1.6% for a level of 8.0 mg / dL.

**Total Bilirubin:** Regression equation  $y = 1.014x - 0.010$  and correlation coefficient  $r = 0.9901$ . Using this equation the total systematic error estimated was 1.20% for a level of 5.0 mg / dL and 1.32% for a level of 13.0 mg / dL.

**Precision:** It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, obtaining:

**Direct Bilirubin**

Samples (mg/dL)	Repetitions	Within-run precision		Total Precision	
		SD (mg/dL)	%CV	SD (mg/dL)	%CV
0.767	80	0.032	4.2	0.039	5.1
3.455	80	0.072	2.1	0.181	5.2
12.218	80	0.067	0.5	0.172	1.4

**Total Bilirubin**

Samples (mg/dL)	Repetitions	Within-run Precision		Total Precision	
		SD (mg/dL)	%CV	SD (mg/dL)	%CV
1.0471	80	0.0135	1.3	0.0135	1.3
5.4074	80	0.0519	1.0	0.0623	1.2
13.4034	80	0.0463	0.3	0.1097	0.8

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

**RESIDUAL RISKS, WARNINGS AND PRECAUTIONS**

- Use PPE according to Good Clinical Laboratory Practice
- Do not mix reagents from different lots.
- Do not exchange reagent bottles to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the reagents out of the specified storage conditions when they are not in use.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

**REFERENCE RANGES**

newborn	Total Bilirubin ( ng/dL )	
Age	Pre mature	Normal
0-1 day	1.0 to 8.0	2.0 to 6.0
1-2 day	6.0 to 12.0	6.0 to 10.0
3-5 day	10.0 to 14.0	4.0 to 8.0

Adults		mg/dL
Direct Bilirubin	0.0 to 0.2	
Total Bilirubin	0.0 to 2.0	

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

**MATERIALS REQUIRED TO CARRY OUT THE TEST**

- Spectrophotometer or photometer for reading at 550 nm (530-570).
- Glass and / or automatic pipettes.
- Cock or Stopwatch.
- Test tubes.

**ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL**

• Dispose, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at [www.bioteccalindbr.com](http://www.bioteccalindbr.com) or by telephone (35) - 3214-4646.

• Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGSS).

**QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE**

Before being approved for use Biotecca reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiry date stated in the label of the external packaging, since it is stored and transported in the specified conditions. Adaptations of semi-automatic protocols with preparation of Working Reagents, data concerning the Quality Control of this product (lot printed on the labels of the reagent bottles) or any doubt in the use of this kit, contact the Scientific Advisory Board of Biotecca Ltda, by phone +55 35 3214 4646 or by email [sac@bioteccaltda.com.br](mailto:sac@bioteccaltda.com.br).

**AUTOMATION**

This procedure is automated in most analyzers. The applications are available at [www.bioteccaltda.ind.br](http://www.bioteccaltda.ind.br)

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