CKMB

CKMB | CKMB

Ref. 11.003.00

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ENGLISH

INTENDED USE



Kit intended to determination of Creatine Kinase MB (subunit B) enzymatic activity in serum. Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

This method is based on the determination of the Creatine Kinase (CK) activity in presence of an antibody specific to the monomer M of CK. This antibody inhibits the isoenzyme CK-MM and half of the CK-MB activity of CK-MB, without affecting the monomer B Isoenzyme CK-MB and CK-BB. As the CK-BB isoenzyme is not found in the blood normally, determining the amount of the monomer B is almost specific for CK-MB form. The catalytic concentration of CK-B matching the half of the CKMB activity is measured by the coupling reactions of Hexokinase and Glucose-6-phosphate Dehydrogenase (G6P-DH), from the speed of the formation of NADPH at 340 nm

SAMPLE – PREPARATION AND STABILITY

Sample Type: Serum

Collection, handling and preparation: Make the collection of the sample in accordance with Good Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: The creatine kinase MB in serum is stable for 7 days if stored at a temperature of 4 to 8 °C, 8 hours at 20 to 25 °C, and 4 weeks at -20 °C. Do not use hemolyzed samples

PRODUCT DESCRIPTION

R 1	imidazole buffer> 40 mM, D-glucose> 10mM, magnesium acetate> 10 mM, NADP> 1.5 mM, N-Acetyl Cysteine> 10 mM, Monoclonal Antibody Anti-CKMM> 2ml / L, Hexokinase> 3,000 U / L, activators, preservative. 
R 2	Creatine phosphate> 50 mM, ADP> 2mM, Di-adenosine-5-pentafosphate> 10 uM, AMP> 10 mM, glucose 6-phosphate dehydrogenase> 2,500 U / L, stabilizers, activators, preservatives. 
CONTROL	CK-MB Control Solution. Phosphate Buffer> 6 mM, Glycerol> 10% v/v, animal tissue extract> 40U / L, preservative. 

The determination of CK MB is traceable to reference material ERM AD455 / IFCC.

STABILITY IN USE

- The stability of the product (R1 and R2) in use is 12 months as long as followed by the recommended storage conditions (2 to 8 °C).
- The stability of working reagent is 14 days, as long as followed by the conditions of preparation and recommended storage (2 to 8 °C).
- The CONTROL stability in use is 30 days, as long as followed by the conditions of preparation and recommended storage (2 to 8 °C).
- Reagents must remain outside the specified temperature only the time required for testing

TECHNICAL PROCEDURE

A) REAGENT PREPARATION

Work Reagent (WR)

Mix in proportion: 4 parts of R1 + 1 part of R2.

Homogenize it gently. The Work Reagent is stable 14 days at 2 to 8 °C.

CKMB CONTROL

The CK-MB control is liquid and ready to use. After opened it is stable up to 30 days at 2 to 8 °C.

B) OPERATING RANGE

The product operating range is from 7.16 U / L to 600 U / L. For higher values, dilute the sample with NaCl 150 mM (0.9%), repeat the assay and multiply the result by the dilution factor.

QUALITY CONTROL

Use of controls should be a routine practice in the lab. It is suggested to use a control within the reference range or decision level and another control with value in another range of clinical significance. For Laboratory Quality Internal Control is recommended the use of CK-MB control provided in the kit.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pre-heat the reagent for 3 minutes at 37 °C.

2. Pipette

Sample	40 µL
Working Reagent	1,0 mL

3. Mix and transfer immediately to a cuvette. Put it in the cuvette holder thermostated. Trigger the timer.

4- After 3 minutes, record the initial absorbance (A0) and make new readings after exactly 1, 2 and 3 minutes (A1, A2 and A3) respectively.

B) CALCULATIONS

Using the absorbance readings calculate the mean variation of absorbance per minute (ΔA / min):

$$\Delta A/\text{min} = \frac{(A1 - A0) + (A2 - A1) + (A3 - A2)}{3}$$

CK-MB activity in the sample is calculated by multiplying ΔA/min, using the following factor

$$\text{CK-MB (U/L)} = \Delta A/\text{min} \times 8254$$

Exemplo:
 A0 = 0,934 A1 = 0,965
 A2 = 0,996 A3 = 1,027

$$\Delta A/\text{min} = \frac{(0,965 - 0,934) + (0,996 - 0,965) + (1,027 - 0,996)}{3}$$

$\Delta A/\text{min} = 0,031$ CKMB (U/L) = $0,031 \times 8254 = 256$ U/L

C) INTERPRETATION

The given result corresponds to the activity of CK-MB isoenzyme in the sample. According to the Guidelines of the Brazilian Society of Cardiology, in the course of an acute myocardial infarction occurs the elevation of biochemical markers corresponding to the damage. These markers are a set of macromolecules released into the bloodstream when myocardial cells are damaged, occurring loss of membrane integrity of such cells and the dissemination of their enzymes in the interstitial space, through the capillaries and lymph vessels. The creatine kinase MB (CK-MB) is a marker traditionally used being useful in confirming the diagnosis and prognosis aid, because there is a direct association between elevated serum markers and the risk of cardiac events in the short and medium term.

INTERFERING AND LIMITATIONS

Anticoagulants: All anticoagulants may interfere in the reaction.

Hemolyzed, Jaundice and Lipemic Sera: Hemoglobin > 25 mg/dL, Bilirubin > 35 mg/dL, Triglyceride > 500 mg/dL.

PERFORMANCE CHARACTERISTICS

Sensitivity: Detection limit: 4.37 U / L / quantification limit: 7.16 U / L.

Analytical Specificity: The product determines CKMB specifically in the presence of other interfering substances in the sample, up to the concentrations stated above

Accuracy: The method was compared with a similar method for the determination of 40 samples in duplicate. The regression equation obtained was $y = 0,993x + 1.6$ and the correlation coefficient $r = 0,9988$. Using this equation the total systematic error estimated is 4.63% to a level of 300 U / L and 0.9% to a level of 100 U / L.

Precision: Determined using samples in 03 levels of decision, with two runs in duplicate per day during 20 days, being obtained.

Samples (U/L)	Repetitions	Within-Run Precision		Total Precision	
		SD	%CV	SD	%CV
26,23	80	0,51	2,0	0,51	2,0
152,84	80	3,17	2,1	3,45	2,3
448,89	80	4,26	0,9	5,1	1,1

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with Good Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any sign in disagreement with the one specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions, when they are not in use.
- Use glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The level of the water bath must be greater than the test tubes containing the reaction.
- The level of water in the thermostatic bath must be over the level of mix reaction within the assay tubes.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

REFERENCE RANGES

Temperature	37 °C
ADULTS	< 24 U/L

These values are intended for orientation only. It's recommended that each lab establishes its own reference range.

Conversion for International System of Units (SI):

CK-MB (U/L) X 0,017 = CK-MB (μKat/L).

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer with a thermostated cuvette able to read at 340 nm.
- Thermostatic water bath at 37°C.
- Glass Pipettes and/or automatic.
- Clock or chronometer.
- Assay tubes

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotecnica.ind.br or calling for 55-35-3214-4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sae@biotecnicaltda.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.biotecnicaltda.ind.br

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- SOCIEDADE BRASILEIRA DE CARDIOLOGIA. Diretrizes Sobre Angina Instável e Infarto Agudo do Miocárdio Diretrizes de Prevenção da Aterosclerose. **Arq. Bras. Cardiol.** v.102, supl. 1, 2014.
- SCHUMANN, G.; et al. IFCC Primary Reference Procedures for the Measurement of Catalytic Activity Concentrations of Enzymes at 37°C – Part 2. Reference Procedure for the Measurement of Catalytic Concentrations of Creatine Kinase. **Clin Chem Lab Med** , v. 40(6), p. 635-642, 2002.
- MYOCARDIAL INFARCTION REDEFINED – A Consensus Document of the Joint European Society of Cardiology/ American College of Cardiology/ American College of Cardiology Committee for the Redefinition of Myocardial Infarction. *Eur Heart J.* v. 21, p. 1502-1513, 2000.
- BURTIS, C. A.; ASHWOOD, E. R. Tietz textbook of clinical chemistry, 4 ed. Philadelphia: W.B.Saunders, 1994, 836p.
- YOUNG, D. S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.

TABELA DE SÍMBOLOS INTERNACIONAIS / TABLE OF INTERNATIONAL SYMBOLS / TABLA DE SÍMBOLOS INTERNACIONALES

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