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
**ENGLISH**  
**INTERFERENCE**  
Kit intended to determine Total Creatine Kinase activity in serum and plasmas. Diagnostic use only.

**STORAGE AND HANDLING**

- Store at 2 to 8°C
- Protect from light.

- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired

**WORKING PRINCIPLE**  
The enzyme Creatine Kinase (CK) catalyses the phosphorylation of ADP to ATP by the phosphate from creatine phosphate. The activity is then determined by coupling the reactions of Hexokinase and Gucose-6-phosphate Dehydrogenase detected by the velocity of the formation of NADPH at 340 nm

**SAMPLE - PREPARATION AND STABILITY**  
**Sample Type:** Serum  
**Collection, handling and preparation:** Make the collection of the sample in accordance with Good Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.  
**Preservation:** The creatine kinase in serum is stable for 30 days if stored at a temperature of 4 to 8°C or -20°C and 4 hours at 20°C to 25°C. Keep samples protected from light.

**PRODUCT DESCRIPTION**

<b>R 1</b>	Indicad e buffer >50 mmol/L, NADP > 1 mmol/L, D-glicose > 20 mmol/L, N-Acetyl-cisteine > 20 mmol/L, Hexokinase > 2000 U/L, activators, preservative, Creatine phosphate > 100 mmol/L, AMP > 5 mmol/L, D-adenosina 5-pentafosfato > 20 µmol/L, ADP > 5 mmol/L, glicose 6-phosphate dehidrogenase > 4000 U/L, activators, preservative.	<b>X</b>
<b>R 2</b>		<b>X</b>

The determination of CK is traceable to reference material IFCC

**STABILITY IN USE**

- The stability of the product (R1 and R2) in use is 15 months as long as followed by the recommended storage conditions (2 to 8°C).
- The stability of working reagent is 14 days, as long as followed by the conditions of preparation and recommended storage (2 to 8°C).
- Reagents must remain outside the specified temperature only the time required for testing

**TECHNICAL PROCEDURE**  
**A) REAGENT PREPARATION**  
**Work Reagent (WR)**  
Mix in proportion of 4 parts of R1 + 1 part of R2.  
Homogenize gently. The Work Reagent is stable 14 days at 2 to 8°C  
**B) OPERATING RANGE**  
The product operating range is from 8,35 U/L to 1000 U/L. For higher values, dilute the sample with NaCl 150 mM (0,9%), repeat the assay and multiply the result by the dilution factor.

**QUALITY CONTROL**  
Use of controls should be a routine practice in the lab. It is suggested to use a control within the reference range or decision level and another control with value in another range of clinical significance. For the internal quality control of the lab it is indicated the use of the calibrator serum and control sera below

Soro Calibrador - Autocal H REF 13.002.00  
Soro Controle Normal - Quanti nor m REF 13.003.00  
Soro Controle Patológico - Quanti at REF 13.004.00

**TEST PROCEDURE, CALCULATION AND INTERPRETATION**  
**A) TEST PROCEDURE**  
1. Heat the work reagent for 3 minutes at 37°C  
2. Zero the reading equipment with purified water filter at 340 nm  
3. Pipette in the assay tube:

	Volume
Work Reagent	1.0 mL
Sample	40 µL

4. Homogenize and insert it immediately in the thermostated cuvette at 37°C  
5. After 2 minutes note the initial absorbance (A0) and read again after exactly 1, 2 and 3 minutes (A1, A2 and A3, respectively).

**B) CALCULATION**  
Using the measured absorbances, calculate the mean variation of absorbance per minute (ΔA/min):  
(ΔA/min) =  $\frac{(A1 - A0) + (A2 - A1) + (A3 - A2)}{3}$

The activity of the CK is then calculated by multiplying the ΔA/min by the factor below

$$CK (U/L) = \Delta A / \text{min} \times 4127$$

Exemplo: A0 = 1,165 A1 = 1,178 A2 = 1,190 A3 = 1,199  
ΔA / min =  $\frac{(1,178 - 1,165) + (1,190 - 1,178) + (1,199 - 1,190)}{3}$

$\Delta A / \text{mín} = 0.011$   
 $CK (U/L) = 0.011 \times 4127 = 45.4 U/L$   
 Note: To calculate the Factor for automatic analyzers use the absorption coefficient of NADPH = 6.3

**Q INTERPRETACION**  
 CK is a diagnostic enzyme that occurs in four different forms: a mitochondrial isoenzyme and the cytosolic isoenzymes CK-MM (skel. muscle), CK-BB (brain) and CK-MB (myocardial). The determination of CK activity and its isoenzymes are used in the diagnosis and monitoring of myocardial infarction and myopathies such as Duchenne muscular dystrophy. In neurological muscular diseases there is no change of serum enzyme activity. It is very high in malignant hyperthermia.

**I INTERFERENCIAS Y LIMITACIONES**  
**Anti coagulantes:** All anticoagulants may interfere in the reaction.  
**Hemolized, Jaundiced and Lipemic Sera:** Hemoglobin > 200 mg/dL, Bilirubin > 30 mg/dL, Triglycerides > 2000 mg/dL.  
**Drugs:** see recommended reference (Young, 2000).

**PERFORMANCE CHARACTERISTICS**  
**Sensitivity:** Detecti on limit: 4.43 U/L / quantifi cation limit: 8.35 U/L  
**Analytical Specificity:** The product determines CK specifically in the presence of other interfering substances in the sample, up to the concentrations stated above.  
**Accuracy:** The method was compared with a similar method for the determination of 40 samples in duplicate. The regression equation obtained was  $y = 1.020x - 3.48$  and the correlation coefficient  $r = 0.9951$ . Using this equation the total systematic error estimate is -6.71% to a level of 40U/L and 0.06% to a level of 180 U/L.  
**Predisi on:** Determined using samples in 0.3 levels of decision, with two runs in duplicate per day during 20 days, being obtained:

Samples (U/L)	Repetitions	Within-Run Prediction		Total Prediction	
		SD (U/L)	%CV	SD (U/L)	%CV
176,178	80	1,502	0,9	2,302	1,3
269,535	80	5,864	2,2	7,050	2,6
74,394	80	1,756	2,4	2,279	3,1

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

**RESIDUAL RISK, WARNINGS AND PRECAUTIONS**

- Use protective equipment in accordance with Good Laboratory Practice
- Do not mix reagents from different lots.
- Do not exchange the caps from different reagents in order to avoid cross contamination
- Do not use the reagent if it displays any significant disagreement with the one specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions, when they are not in use.
- Use glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The level of the water bath must be greater than the test tubes containing the reaction.
- The level of water in the thermostatic bath must be over the level of mix reaction within the assay tubes.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

**REFERENCE RANGES**

	37 °C
ME N	< 170 U/L
WO ME N	< 190 U/L

These values are intended for orientation only. It's recommended that each lab establishes its own reference range.  
 Conversion for International System of Units (SI):  
 $CK (U/L) \times 0.017 = CK (\mu kat)$

**ADDITIONAL EQUIPMENT**

- Spectrophotometer or photometer with a thermostated cuvette able to read at 340 nm
- Thermostatic water bath at 37°C
- Glass pipettes and/or automatic
- Clock or chronometer.
- Assay tubes

**WARNINGS AND PRECAUTIONS**

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infective material.

- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at [www.bioteccaltda.com.br](http://www.bioteccaltda.com.br) or calling for 55-35-3214-4646

**QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE**

Before being approved for use, the biochemical reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiry date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the label of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: [sac@bioteccaltda.com.br](mailto:sac@bioteccaltda.com.br)

**AUTOMATION**

This product is compatible to the most types of biochemical automatic analyzers. The applications are available at [www.bioteccaltda.com.br](http://www.bioteccaltda.com.br)

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- THOMAS, L.; et al. Consensus of DGKL and VDGHG for interreference intervals on enzymes in serum *J Lab Med*, v. 29, p. 301-308, 2005.
- SCHUMANN, G.; et al. IFCC Primary Reference Procedures for the Measurement of Catalytic Activity Concentrations of Enzymes at 37°C. Part 2. Reference Procedure for the Measurement of Catalytic Concentration of Creatine Kinase. *Clin Chem Lab. Med.* v. 40 p. 635-642, 2002.
- BURTIS, C. A.; ASHWOOD, E. R.; BRUNS, D. E. *Tietz Fundamentos de Química Clínica*, Saunders Elsevier, 6 ed. 2008
- YOUNG, D.S. *Effects of drugs on clinical laboratory tests - vol. 2*, 5 ed. Washington DC: AACCPress, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin Chem* v 27 p. 493-501, 1981.

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