# CÁLCULO

Calculo | Cálculo

Ref. 12.001.00

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**ENGLISH**

**INTENDED USE**  
Kit intended to determine the concentration of calcium in serum, plasma and urine. Diagnostic only.

**STORAGE AND HANDLING**  
• Store at 15 to 30 °C  
• Protect from light.  
• Stable until the kit expiration date that is printed on the label.  
• Do not use reagents whose date has expired.

**WORKING PRINCIPLE**

Calcium ions with orthocresolphthaléin complexone, in alkaline medium, a colored complex that can be measured at 550-590 nm. The intensity of the formed color is proportional to the amount of calcium present in the sample.

**SAMPLE - PREPARATION AND STABILITY**

**Sample and Type:** Serum, plasma (heparin) and urine.

**Collection, handling and preparation:** Make the collection of the sample in accordance with Good Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.  
**Preservation:** The calcium in serum and plasma (heparin) is stable for 14 days if stored at a temperature of 4 to 8 °C and 6 months if stored at a temperature of -20 °C. Uri na: Use sample collected within 24 hours in a bottle containing 20 to 30 mL of 6 m/L HCl solution to prevent the precipitation of calcium salts. Mix well, measure the volume and proceed with the test.

**PRODUCT DESCRIPTION**

<b>R 1</b>	2-Anilino methyl propanol > 1% v/v; Monoethanolamine > 1% v/v and stabilizers.		
<b>R 2</b>	Hydrochloric acid > 0,1% v/v, Orthocresolphthaléin complexone > 0,001% w/v, 8-hydroxyquinoline > 1 mmol/L and preservative.		
<b>STD</b>	Calcium chloride in concentration equivalent to 10 mg/dL and preservative. Traceable to reference material NIST 3109a.		

**STABILITY IN USE**

- The stability of the product (R1, R2 and STD) in use is 24 months as long as followed by the recommended storage conditions (15 to 30 °C).

**TECHNICAL PREPARATION**

**A) REAGENT PREPARATION**  
**Working Reagent:** Mix in proportion: 1 part of R1 + 1 part of R2. Homogenize gently. Stable for 24 hours at 15 to 30 °C.  
**NOTE:** The absorbance of the working reagent should be less than 0,200 at 570 nm in equipment cleared with purified water.

**STD Ready for use**  
**B) OPERATING INTERVAL**  
The operating range of the product is from 0,91 ng/dL to 18,00 mg/dL. For higher values, dilute the sample with 150 mM NaCl (0,9%), carry out a new dosage and multiply the result obtained by the dilution factor.

**QUALITY CONTROL**

Use of controls should be a routine practice in the lab. It is suggested to use a control within the reference range or decision level and another control with value in another range of clinical significance. For the internal quality control of the laboratory indicated the use of the calibrator serum control should be low.  
Calibrator Serum- Autocal H 13.002.00  
Normal Control Serum- Quantinorm 13.003.00  
Pathologic Control Serum- Quantinorm 13.004.00

**TEST PROCEDURE, CALCULATION AND INTERPRETATION**

**A) TEST PROCEDURE**

1. Pipette in the assay tubes:

	Blank	Standard	Sample
STD	-	10 µL	-
Sample	-	-	10 µL
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate 2 minutes at 15 to 30 °C
3. Measure the absorbance of Standard (Ast) and Sample (Asa) against the blank at 570 nm (550-590 nm). The final reaction is stable for 30 min.

**B) CALCULATIONS**

**Calculation of serum**  
Cálculo (ng/dL) =  $\frac{\text{Sample Absorbance}}{\text{Standard Absorbance}} \times \text{STD Concentration}$

**Exemplo**  
Standard Concentration = 10 ng/dL  
Sample Absorbance = 0,220  
Standard Absorbance = 0,250  
Cálculo (ng/dL) =  $\frac{0,220}{0,250} \times 10 = 8,8$  mg/dL

**With Calibration Factor (CF):**

CF =  $\frac{\text{STD Concentration (ng/dL)}}{\text{STD Absorbance}}$   
Cálculo (ng/dL) = Sample Absorbance x CF

**Exemplo**

Calibration Factor =  $\frac{10}{0,250} = 40$   
Cálculo (ng/dL) = 0,220 x 40 = 8,8 mg/dL

**Uri na 24 hours:**

Calcium (mg / dL) =  $\frac{\text{Sample Absorbance} \times \text{Standard Concentration (ng / dL)}}{\text{Standard Absorbance}}$   
 Calcium (ng/24h) =  $\frac{\text{urine calcium (ng/dL)} \times \text{urinary volume of 24 hr (ml)}}{100}$

**Q I NTERPRETATI ON**

Calcium is the most predominant cation in the body. The skeleton contains 99% of body calcium with an unknown structure and composition roughly that of hydroxyapatite. The remainder of the calcium is distributed in tissues and extracellular fluid where it performs its specific functions in vital processes, among them: blood coagulation, neuromuscular conduction, cardiac and skeletal muscle excitability, enzymatic activation and the preservation of cell membrane permeability and integrity. Serum calcium levels are controlled by parathyroid hormone (PTH), calcitonin and vitamin D. An imbalance in some of these modulators causes changes in serum and body calcium levels. Serum increases in PTH or vitamin D are generally associated with hypercalcemia. Increased levels of serum calcium can be observed in multiple myeloma and other neoplastic diseases. Hypocalcemia can be observed for example in hypoparathyroidism, nephrosis and pancreatitis.

**I NTERFERENCI A S Y I NTERFERENCI AS**

**Hemolysis, Jaundice and Lipemic Serum:** Hemoglobin > 500 mg/dL, Bilirubin > 35 mg/dL, Triglyceride > 1000 mg/dL interfere in the dosage.

**Drugs:** see recommended reference (Young, 2000).

**PERFORMANCE CHARACTERISTICS**

**Sensitivity:** Detection limit: 0.58 mg/dL / quantification limit: 0.91 mg/dL  
**Analytical Specificity:** The product determines calcium specifically in the presence of other interfering substances in the sample up to the concentrations stated above.  
**Accuracy:** The method was compared with a similar method for the determination of 40 samples in duplicate. The regression equation obtained was  $y = 0.990x + 0.090$  and the correlation coefficient  $r = 0.9984$ . Using this equation on the total systematic error estimate is 0.80% to a level of 5.0 mg/dL and -0.25% to a level of 12.0 mg/dL.  
**Precision:** Determined using samples in 03 levels of decision, with two runs in duplicate per day during 20 days, being obtained:

Samples (mg/dL)	Repetitions	Within-Run Precision		Total Precision	
		SD	%CV	SD	%CV
10,2688	80	0,0564	0,5	0,1332	1,3
13,2685	80	0,0229	0,2	0,1021	2,1
16,9265	80	0,0118	0,1	0,1542	0,9

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

**RESIDUAL RISK WARNINGS AND PRECAUTIONS**

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator, and reagent.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meet them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are key factors for the stability of the reagents and obtaining correct results.

**REFERENCE RANGES**

	mg/dL	mmol/L
Serum and Plasma	8,6 - 10,3	2,15 - 2,57
	ng/24 hours	mmol/24 hours
Urine	Up to 300	Up to 7,49

These values are intended for orientation only. It is recommended that each lab establishes its own reference range. Conversion for International System of Units (SI): mmol/L Calcium (mg/dL) X 0.25 = Calcium (mmol/L)

**ADDITIONAL EQUIPMENT**

- Spectrophotometer or photometer able to read at 570 nm (550-590).
- Pipettes and micro pipettes.
- Clock or chronometer.
- Assay tubes.

**WARNINGS AND PRECAUTIONS**

- Discard the reaction surplus according to Good Laboratory Practice in a proper place for potentially ineffective material.
- The information on Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at [www.bioteccaltda.com](http://www.bioteccaltda.com) or calling for 55-35-3214-4646

**QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICES**

Before being approved for use, Biotécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiry date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the label)

of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: [sac@bioteccaltda.com.br](mailto:sac@bioteccaltda.com.br)

**AUTOMATION**

This product is compatible to the most types of biochemical automatic analysers. The applications are available at [www.bioteccaltda.com.br](http://www.bioteccaltda.com.br)

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