IRON FERROZ INE BT 12.005.00

UBC = Latent Iron Binding Capacity
TIBC = Total Iron Binding Capacity
TSI = Transferrin Saturation Index

UBC (µg/dL) = 500 - [Standard Absorbance X 500]
TIBC (µg/dL) = UBC + serum iron
TSI (%) = Serum Iron X 100 / TIBC

Exa mple
A1 sample = Q 042 Standard Absorbance = Q 480
A2 sample = Q 391 Serum iron: 110 µg/dL

UBC (µg/dL) = 500 - [(Q 391 - Q 042) / Q 480 X 500] = 136 µg/dL
TIBC (µg/dL) = 110 + 136 = 246 µg/dL

TSI (%) = 110 X 100 = 44.7% / 246

SENSITIVITY AND LINEARITY
Methodological sensitivity (UBC)
Calculated: 1,075 µg/D
Linearity: 450 µg/dL
For values higher, dilute the sample with 0.9% NaD 150 mM solution, repeat the assay and multiply the obtained result by the dilution factor.

TECHNICAL LIMITATIONS
Anti coagulants
All the anti coagulant interfere in the reaction
Jaundice and Lipemic Sera
Total Bilirubin > 10 mg/dL
Lipemic > 900 mg/dL

QUALITY CONTROL
Any clinical laboratory must keep an internal quality control program which defines the aids, procedures and criteria for the tolerance limits, corrective actions and registration of the activities. Also it must be kept a defined system for verification of the analytical variability that occurs in any measuring system.
The use of controls to evaluate the imprecision of the analysis must be a routine practice in the lab. It's suggested to use a control within the reference range and other control within the critical significance range. The application of the Multiple Rules of Westgard for evaluation of the control status is recommended.
The lab must participate of external quality control programs.
For the internal quality control of the lab it's indicated the use of the calibrator serum and control sera below
Normal Control Serum - Quanti nor 13.003.00
Pathological Control Serum - Quanti Art REF 13.004.00

REFERENCE RANGES
Table with 2 columns: Parameter (UBC, TIBC, TSI) and Range (140 to 280 µg/dL, 250 to 410 µg/dL, 20 to 50 %)

These values are intended for orientation only. It's recommended that each lab establishes own reference range.

Conversion for International System of Units (SI): mmol/L
TIBC (µg/dL) X 0.179 = mmol/L

PERFORMANCE CHARACTERISTICS
Intra Assay
The realization of 20 determinations of the same sample at the same day showed a Coefficient of Variation of 4.21%
Inter Assay
The realization of 10 determinations of the same sample at different days showed a Coefficient of Variation of 4.46%
Analytical Specificity
A comparison with a reference method showed a correlation coefficient (r) of 0.989 obtained from a non linear correlation. The result equation of the linear regression is Y = 1,0041 X + 7,8468

CLINICAL SIGNIFICANCE
The plasma iron transport is carried out by a specific beta-globulin the Transferrin which may be determined by the iron quantity it can bind. This quantity is called the Total Iron Binding Capacity (TIBC).
For the complete diagnostic of accumulated iron or deficiency of iron in the patient, it's important to perform a complete haematological profile, iron determination, TIBC and Transferrin Saturation Index (TSI).
The frequent reasons for the increasing TIBC are: anaemia by iron deficiency, pregnancy, use of contraceptives, hypoxia events, etc. The TIBC is decreased in chronic inflammatory diseases, malignant diseases, nephropathies, iron intoxication, talasemia major. The TSI is increased in iron intoxication, nephropathies, talasemia major, pyridoxine deficiency, haemolytic and sideroblastic anaemias. The TSI is decreased in iron deficiency, malignant diseases, chronic infections, pregnancy and menstruation periods.

OBSERVATIONS
1. The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be

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validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
2. The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

QUALITY ASSURANCE/CUSTOMER TECHNICAL SERVICE
Before being approved for use, Biotécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiry date stated in the label of the external packaging, since it is stored and transported in the specified conditions.
The quality control data concerning this product or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending us an e-mail: sac@biocnica.com.br

AUTOMATION
This product is compatible to the most types of biotechnical automatic analysers. The applications are available at: www.biocnica.com.br

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- GOODWIN J.; MURPHY, B.; GUILLETTE M Direct measurement of serum iron and binding capacity. Clin Chem v.12, p 47-57, 1966.
STOCKEY, L. L. Ferronine - a new spectrophotometric reagent for iron. Anal. Chem v.42, p 779-781, 1970.
ITANQ - M CAP Comprehensive chemistry: serum iron survey. Am J. Clin Pathol. v.70, p516-522, 1978.
YOUNG D.S Effects of drugs on clinical laboratory tests - v.d. 2 5ed. Washington DC: AACC Press, 2000.
WESTGARD J. O et al. A multi-rule shewhart chart quality control in clinical chemistry. Clin Chem v.27 p 493-501, 1981.

TABELA DE SÍMBOLOS INTERNACIONAL / TABLE OF INTERNATIONAL SYMBOLS / TABLA DE SÍMBOLOS INTERNACIONALES
Table with 4 columns: Symbol (i, REF, LOT, IVD, R <N>), Description (Consult instructions, Code, Batch code, Reagent, etc.), Symbol (Σ, pipette, flask, X, STD), and Note (Control sufficient, etc.)