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ENGLISH

INTENDED USE

Kit intended to detect concentration of cholesterol in serum Diagnostics only.


STORAGE AND HANDLING

- Store at 2 to 8°C
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

Cholesterol esters from sample are hydrolyzed by the enzyme cholesterol esterase producing free cholesterol. The enzyme cholesterol oxidase, in presence of oxygen,

catalyses the oxidation of the free cholesterol producing hydrogen peroxide. The enzyme peroxidase catalyses the oxidation of phenid by the hydrogen peroxide in presence of 4-Aminoantipirine, producing a pink compound (4-(p-benzoquinononil)fenazon) which is detected in the spectrophotometer at 505 nm

2 H₂O + 4 Aminoantipirine + Phenid $\xrightarrow{\text{peroxidase}}$ 4-(p-benzoquinononil)fenazon + 4H₂O

SAMPLE - PREPARATION AND STABILITY

Sample Type: Serum

Collection, handling and preparation: Make the collection of the sample in accordance with Good Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: The Cholesterol is stable in serum for 7 days from 4 to 8°C and 3 months from -20°C

PRODUCT DESCRIPTION

R 1 

STD 

pH buffer \geq 20 mmol/L; 4 aminoantipirine \geq 0,1 mmol/L; phenid \geq 1,0 mmol/L; peróxido \geq 1000 U/L; chd ester oxidase \geq 200 U/L; chd ester esterase \geq 200 U/L; lipoproteína lipase \geq 200 U/L; activators; detergent; stabilizers and preservative; Butylated chd \geq 5% v/v, detergent, stabilizers and Cholesterol concentration of 200 mg/dL. Traceable to reference material NST 1951b

STABILITY IN USE

- The stability of the reagents R1 and STD in use is 18 months as long as followed by the recommended storage conditions (2 to 8°C).
- Reagents must remain outside the specified temperature only the time required for testing

TECHNICAL PROCEDURE

A) REAGENT PREPARATION
Reagent 1 e Standard (R1 / STD)

Ready-to-use

Note: Due to Standard validity, the bottle must be decanted after its use

B) OPERATING RANGE

The product operating range is from 4,60 ng/dL to 800,00 mg/dL. For higher values, dilute the sample with NaCl 150 mM (0,9%), repeat the assay and multiply the result by the dilution factor.

QUALITY CONTROL

Use of controls should be a routine practice in the lab. It is suggested to use a control within the reference range or decision level and another control with value in another range of clinical significance. For the internal quality control of the laboratory used the use of the calibrator serum control sera below

Calibrador Serum- Autocal REF 13.002.00

Normol Control Serum- Quantinorm REF 13.003.00

Padrão de Colesterol Control Serum- Quantialt REF 13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

- Pipette in the assay tubes:

	Blank	Standard	Sample
STD	-	10 µL	-
Sample	-	-	10 µL
R1	1,0 mL	1,0 mL	1,0 mL

- Homogenize and incubate 10 minutes at 37°C
- Measure the absorbance of Standard (Ast) and Sample (Asa) against the Blank at 505 nm. The final reaction is stable for 25 minutes.

B) CALCULATION
Calculation of Serum

Colesterol (ng/dL) = Sample Absorbance x STD value (stated on the label) / Standard Absorbance

With Calibration Factor (CF):
CF = STD Concentration (ng/dL) / STD Absorbance

Colesterol (ng/dL) = Sample Absorbance x CF

C) INTERPRETATION

Cholesterol is a steroid with a secondary hydroxyl group on the C3 position and is an essential component of the membrane of all cells. Approximately 75% of it is synthesized endogenously in the liver and intestine, the remainder is obtained through diet. Most peripheral cells depend on the exogenous distribution of cholesterol by the lipoproteins. Its elevation may be associated with atherosclerotic risk or other metabolic disorders of lipids and lipoproteins.

INTERFERING AND LIMITATIONS

Anticoagulants: All anticoagulants may interfere in the reaction

He mol yzad, Jaundic e and Li penic Será: Hemógl obin > 200 ng/dL, Bilirrubin > 15 mg/dL, Triglicerid e > 1250 ng/dL.

Dr ugs: see recommended reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Detecti on limit: 1,74 ng/dL / quantificati on limit: 4,60 ng/dL.

Analytical Specificity: The product determines Cholesterol specifically up to the concentrations of interfering substances in the sample, stated above.

Accuracy: The method was compared with a similar method for the determination of 40 samples in duplicate. The regression equation obtained was $y = 0,992x - 4,48$ and the correlation coefficient $r = 0,9908$. Using this equation the total systematic error estimate is $-3,04\%$ to a level of 200 ng/dL and $-1,69\%$ to a level of 500 ng/dL.

Precision: Determined using samples in 03 levels of decision, with two runs in duplicate per day during 20 days, being obtained:

Samples (ng/dL)	Repetitions	Within-Run Precision		Total Precision	
		SD (ng/dL)	%CV	SD (ng/dL)	%CV
135,77	80	1,20	0,9	1,36	1,0
295,52	80	2,45	0,8	2,98	1,0
595,42	80	2,62	0,4	2,86	0,5

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with Good Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any sign of disagreement with the one specified in the product MSDs.
- Avoid leaving reagents outside the specified storage conditions, when they are not in use.
- Use glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The level of the water bath must be greater than the test tubes containing the reaction.
- The level of water in the thermostatic bath must be over the level of mix reaction with in the assay tubes.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

REFERENCE RANGES

According to the Brazilian directive Cholesterol levels can be evaluated based on age as follows:

Age	Total Cholesterol Values	
	(ng/dL)	Interpretation
2 to 19 years	< 150	Desirable
	150 – 169	Border Range
	> 170	High
20 years or more	< 200	Desirable
	200 – 239	Border Range
	> 240	High

These values are intended for orientation only. It's recommended that each lab establishes its own reference range.

Conversion for International System of Units (SI): mmol/L

Cholesterol (ng/dL) X 0.026 = mmol/L

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer with a thermostated cuvette able to read at 505 nm (490 – 510 nm).
- Thermostatic water bath at 37°C.
- Glass Pipettes and/or automatic.
- Clock or chronometer.
- Assay tubes

WARNINGS AND PRECAUTIONS

- Discard the reaction surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.bioteccalind.com.br or calling for 55-35-3214-4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use, Biotécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiry date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels

of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@bioteccalind.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.bioteccalind.com.br

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