## CREATI NA NA

Creati ri ne | Creati ri na - Ref. 10.007.00

### RNALIDADE

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#### ENGLISH

#### INTERNEED USE

Kit intended to determine creatinine in serum plasma and urine. In-Vitro Diagnostic use only.

#### STORAGE AND HANDLING

- Store at 15 to 30°C
- Protect from light.
- Expiration date that is printed on the package label.
- Do not use reagents whose date has expired.

#### WORKING PRINCIPLE

The creatinine in the sample reacts with the picrate in alkaline medium to give an orange-red complex. The intensity of the formed color is proportional to the creatinine concentration in the sample.

#### SAMPLE - TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

**Sample Type:** Serum plasma (EDTA and heparin) and urine.

**Collection, handling and preparation:** Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

**Preservation:** Serum and plasma creatinine is stable for 7 days if stored at 20 to 25°C, 7 days if stored at 4 to 8°C and 3 months at -20°C. Urine creatinine is stable for 2 days if stored at 20 to 25°C, 6 days if stored at 4 to 8°C and 6 months if stored at -20°C. Collect urine for 24 hours. Homogenize the urine, measure the volume and separate a 20 mL sample. Centrifuge for 10 minutes at 3000 rpm. Dilute an aliquot of the centrifuged urine at the ratio of 1:25 with purified water and proceed with the assay. Multiply the result obtained by 25. This dilution should result in a measurement within the operating range, to higher values perform new dilution by changing the ratio.

#### PRODUCT DESCRIPTION

<b>R 1</b>	Picric acid ≥ 10 mmol/L	
<b>R 2</b>	Sodium Hydroxide > 200 mmol/L, acetic acid > 10 mmol/L, detergent.	
<b>STD</b>	Stabilizer; preservative; Creati ri ne (the concentration is equivalent to 2 mg/dL. The determination of creati ri ne is traceable to reference material NIST 914a.	

#### STABILITY IN USE

- After opening the product (R1, R2 and STD) in use is stable to the expiration date printed on the label since the recommended storage conditions (15 to 30°C) are followed.
- The stability of the working reagent is 21 days, since the recommended conditions of preparation and storage are followed: temperature of 2 to 8°C in a plastic bottle closed amber and protected from light. The product should remain outside the specified temperature only time enough to run the tests.

#### TREATMENT AND HANDLING OF THE PRODUCT

##### A) REAGENT PREPARATION

##### Working reagent:

Mix in proportion 1 part of R1 + 1 part of R2. Homogenize gently.

STD: Ready to use

##### B) OPERATIONAL INTERVAL

Serum/Plasma: The operating range of the product is from 0,28 mg/dL to 12,00 mg/dL.

For higher values, dilute the sample with 150 mM NaCl (0,9%), carry out a new dosage and multiply the result obtained by the dilution factor.

Urina: The operating range of the product is from 3,15 mg/dL to 321,35 mg/dL corrected by the dilution factor.

#### QUALITY CONTROL

Use of controls should be a routine practice in the lab. It is suggested to use a control within the reference range or decision level and another control with value in another range of clinical significance. For the internal quality control of the lab it's indicated the use of the calibrator serum and control sera below

Calibrador Serum - Autocal H

Nor mal Control Serum - Quant nor m

Pathológico Control e Serum - Quanti dt

<b>REF</b>	13.002.00
	13.003.00
	13.004.00

#### TEST PROCEDURE, CALCULATION AND INTERPRETATION

##### A) TEST PROCEDURE

- 1) Pre-heat the Work Reagent for 3 min. at 37°C
- 2) Zero the equipment at 500nm (490 - 510 nm) with purified water.
- 3) Pipette in the assay tubes:

	<b>Volume</b>
<b>Work Reagent</b>	1,0 mL
<b>Sample or STD</b>	100 µL

4. Homogenize and insert into the thermostated cuvette at 37°C. Set the stopwatch.

5. Record the absorbance at 30 seconds (A1) and 90 seconds (A2) of the sample and the STD

##### B) CALCULATION

##### Serum/Plasma:

Creati ri ne (mg/dL) =  $\frac{(A2 - A1) \text{ sample} \times \text{STD Concentration (mg/dL)}}{(A2 - A1) \text{ STD}}$

##### With Calibration Factor:

Calibration Factor =  $\frac{\text{STD Concentration (mg/dL)}}{(A2 - A1) \text{ of the STD}}$

Creati ri ne (mg/dL) = (A2 - A1 of the Sample) x Calibration Factor

##### Calculation for Uri ne:

Determine the creati ri ne concentration in the urine sample using the same calculation procedure for serum. Multiply the value obtained by 25 (dilution) to obtain the result in mg/dL. This should be used to calculate the creati ri ne value in mg/24h, according to the following equation:

$$\text{Creati ri ne (mg/24 horas)} = \frac{\text{mg/dL} \times \text{urinary volume (in mL)}}{100}$$

##### Creati ri ne in mg/kg/24 hours:

Creati ri ne (mg/kg/24 horas) =  $\frac{\text{Creati ri ne (mg/24 horas)}}{\text{Vt (kg)}}$

Where: U = urine creati ri ne (mg/dL)

S = serum creati ri ne (mg/dL)

VM = minute volume (24 hours urine volume in mL divided by 1440 minutes).

Cálculo de la de Body Surface (SQ) without use of the Nomogram  
 $SC (m^2) = P^{0.725} \times A^{0.725} \times 0.007184$   
 Where: P = weight in Kg  
 A = height in cm

**Clearance**  
 Clearance (ml/ min) =  $\frac{Debug (ml/ min) \times 1.73}{Body Surface}$

**CLINICAL INTERPRETATION**  
 Creatinine is the final product of phosphocreatine decomposition. It is produced endogenously and released into body fluids at a constant rate and its plasma concentration is maintained within narrow limits, predominantly by glomerular filtration. Consequently, both the plasma creatinine concentration and its renal clearance (creatinine clearance) have been used as markers of glomerular filtration rate.

**INTERFERENCES AND LIMITATIONS**  
**He molysis, Jaundice and Lipemia:**  
 Serum / Plas ma: He molgi ob n> 500 ng / dl / Bilirubi n> 10 mg / dl / Tri glyceri des> 500 mg / dl interferences in the dosage.  
 Uri ne: He molgi ob n> 180 mg / dl / Bilirubi n> 6 ng / dl.  
**Medications:** consult recommended literature reference (Young, 2000).

**PERFORMANCE CHARACTERISTICS**  
**Sensitivity:**  
 Serum / Plas ma: Detección límite: 0.11 ng / dl / Límite de cuantificación: 0.28 ng / dl.  
 Uri ne: Límite de detección: 0.018 ng / dl / Límite de cuantificación: 3.15 ng / dl, corrected by dilution factor.  
**Analytical Specificity:** The product specifically determines creatinine in the presence of other interfering substances in the sample up to the concentrations reported above.

**Accuracy:**  
 Serum The method was compared to a similar method by the determination of 40 samples in duplicate. The regression equation is  $y = 0.998x + 0.040$  and correlation coefficient  $r = 0.9958$  were obtained. Using this equation the estimated total systematic error of 3.8% for a level of 1.0 ng / dl and 0.8% for a level of 4.0 ng / dl.  
 Uri ne The method was compared to a similar method by the determination of 20 samples in duplicate. The regression equation is  $y = 1.002x - 2.343$  and correlation coefficient  $r = 0.9995$  were obtained. Using this equation the estimated total methodological error of 0.04% for a level of 1500 ng / 24h and 0.08% for a level of 2000 ng / 24h.

**Preis on:**  
 Serum / Plas ma: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, being obtained.

Sample (ng/dL)	Repetitions	Within Run Precision		Total Precision	
		SD (ng/dL)	%CV	SD (ng/dL)	%CV
0.9620	80	0.0043	0.4	0.0048	0.5
4.0730	80	0.0123	0.3	0.0168	0.4
10.4144	80	0.1246	1.2	0.1320	1.3

% CV: Coefficient of variation expressed as a percentage; SD: Standard deviation

Uri ne: It was determined using samples at 02 decision levels, two runs in duplicate per day for 20 days, being obtained.

**RESIDUAL RISKS / WARNINGS AND PRECAUTIONS**

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different batches.
- Do not exchange reagent bottles and caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator, and reagent.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

**REFERENCE RANGES**

Serum / Plas ma	Men		Women	
	mg/dL	µmol/L	mg/dL	µmol/L
	0.9 - 1.3	80 - 115	0.6 - 1.1	53 - 97
Uri ne	mg/kg/day		µmol/kg/day	
	14 - 26	124 - 230	11 - 20	97 - 177

Clearance (ml/ min)	
Adults	52.1 – 110

These values are for guidance only and it is recommended that each laboratory establishes its own reference range. Conversion to International System Unit (SI): µmol / L  
 Creatinine (mg / dl) x 88.4 = Creatinine (µmol / L)

**MATERIALS REQUIRED TO CARRY OUT THE TEST**

- Spectrophotometer or photometer for reading at 500 nm (490 - 510).
- Water bath thermostated at 37 ° C
- Glass and / or automatic pipettes.
- Cock or Stop watch.
- Test tubes.

**WARNINGS AND PRECAUTIONS**

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.bioteccalind.com or calling for 55-35-3214-4646.

**QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE**  
 Before being approved for use Biotécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiry date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@bioteccalind.com.br

**ADVERTISMENTS**  
 This product is compatible to the most types of biochemical automatic analysers. The applications are available at: www.bioteccalind.com.br

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