ALFA-1 GLICOPROTEÍNA ÁCID DA PLUS

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ENGLISH

INTENDED USE

Kit intended to determine in serum and plasma.

PRINCIPLE OF THE METHOD

The determination is based on a turbidimetric method in which anti-α1-glycoprotein antibody human added glycoprotein, form with the α1-glycoprotein, an insoluble complex providing a turbidity whose intensity is proportional to the concentration of α1-glycoprotein sample, determined in a spectrophotometer at 340 nm.

REAGENT COMPOSITION

R 1	Goat anti bod es anti-human-α1-ad d glycoprotein. Tris buffer pH 8.2 and stabilizer agents.
------------	---

STORAGE AND STABILITY

- After first opening the reagents are stable until the expiration date printed on the label if stored 2-8°C tightly closed.
- Do not use reagents over the expiration date.
- Store at 2 to 8°C
- Do not freeze and protect from light.
- Reagents must remain outside of the specified temperature only the necessary time for the tests.

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer able to read at 340 nm (320-360).
- Thermostatic water bath at 37°C
- Pipettes and micropipettes.
- Clock or chronometer.
- Assay tubes.
- Multicalibrador PLUS Biotécnica - CAT BT 21.011.00

WARNINGS AND PRECAUTIONS

- Perform sample collection in accordance with Good Laboratory Practice.
- The kit is intended only for in vitro diagnostic use.
- The samples must be considered as potentially infective.
- Wear the individual protection equipment accordingly to Good Laboratory Practice.
- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product. The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, also is available on the site www.biotecnica.ind.br or calling for 55-35-3214-4646.
- Do not combine reagents from different batches
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any significant disagreement with the one specified in MSDS of the product.
- Avoid leaving reagents outside of the specified storage conditions, when they are not in use.
- Use glass pipettes and disposable tips for each sample, control, standard/calibrador and reagent.
- The level of water in the thermostatic bath must be over the level of mix reaction within the assay tubes.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

SAMPLE – PREPARATION AND STABILITY

Fresh serum or plasma with heparin or EDTA.
Stable 7 days at 2 to 8°C or 3 months at -20°C.
Discard dilute and/or hemolyzed samples.
Discard the supernatant and conduct the assay.

TECHNICAL PROCEDURE

A) REAGENT PREPARATION

The reagent is ready to use.

Warnings:

The presence of particles, turbidity or absorbance of reagent higher than 0,300 at 340 nm zeroed with purified water indicates deterioration.

B) CALIBRATION CURVE

1. Dilute Multicalibrador PLUS in NaCl 9g/L as follows:

Dilution	1	2	3	4	5	6
Calibrador (µL)	-	10	25	50	75	100
NaCl 9g/L (µL)	100	90	75	50	25	-
Factor	0	0,1	0,25	0,5	0,75	1

Multiply the concentration of α1-glycoprotein Multicalibrador PLUS by the corresponding factor for the concentration of α1-glycoprotein each dilution.

B) PROCEDURE

- Preheat R1 and the spectrophotometer (cuvette holder) at 37°C

2. Reset the equipment with water purified in the filter 340 nm
3. Rpettei na cuvette

	Vol ume
Sampl e/ Cali brator	7 µL
R1	1,0 mL

4. Mix and insert i n t o the cuvette holder ther mostatically.
5. After 4 minutes, measuring absorbance (A) of the sample / cali brator.
6. The concentration of alpha-1 add glycoprotein sample is calculated by i n t e r p e d a t i o n of their absorbance on the cali brati on curve

CALI BRATI ON CALCULATI ON PROCEDURES

- a) Make a calibration curve by plotting the X-axis the concentration of each cali brator diluti on and the Y-axis value of the correspondi ng absorbance (A).
- b) Determini ng the concentrati on of samples by i n t e r p e d a t i n g the absorbance readi ngs (A) obtai ned for the sampl es from the cali brati on curve.

SENSITI VI TY AND LI NEARI TY

Analytical Sensitivity: 1,55 mg/dL
Linearity limit: 250 mg/dL

For higher values, dilute again in the sample with NaCl solution (0,9%), repeat the assay and multiply the obtained result by the dilution factor.

Prozone: The reaction has no prozone effect up to 1000 mg/dL

TECHNICAL LIMITATIONS

Interferences:

- Anti coagulants: They suffer interference from direct oxalate and fluoride.

- Bilirubin > 10 mg/dL; Hemoglobin > 2g/L; Triglycerides > 1100 mg/dL

QUALITY CONTROL

The use of controls to evaluate the imprecision of the analysis must be a routine practice in the lab. It is suggested to use a control within the reference range or decision level and a control with value in another range of clinical significance. For Internal Calibration and Control it is recommended the use of the cali brator serum and control sera below

Multicalibrator PLUS	REF	21.011.00
Multiparameter Control		21.003.00

REFERENCE RANGES

Reference values Alpha-1-acid glycoprotein (mg/dL)	
Adults	50 upto 120

These values are intended for orientation only. It's recommended that each lab establishes its own reference ranges.

PERFORMANCE CHARACTERISTICS

Intra-Assay

The realization of 20 determinations of the same sample at the same day showed a Coefficient of Variation of 1,38%

Inter-Assay

The realization of 10 determinations of the same sample at different days showed a Coefficient of Variation of 2,83%

Analytical Specificity

A comparison with a reference method showed a correlation coefficient (r) of 0,9524 obtained from ambulatory samples. The result equation of the linear regression is $Y = 1,010 X - 1,31$

CLINICAL SIGNIFICANCE

The alpha-1-acid glycoprotein (orosomucoid) is an acute phase protein that has a regulatory effect influencing the inflammatory cascade, thereby protecting the excessive damage on tissues during inflammatory process. It is synthesized by the cells of the hepatic parenchyma, also granulocytes and monocytes contribute significantly to increased plasma levels during sepsis.

The plasma concentration of alpha-1-glycoprotein increases 3-4 times in most conditions associated with inflammation or tissue necrosis, and is one of the most reliable indicators in clinical activity of ulcerative colitis. The levels can also be increased along with haptoglobin and prealbumin by the effect of both endogenous glucocorticoids (Cushing's syndrome) as exogenous.

The reduction in synthesis and plasma levels of alpha-1-acid glycoprotein occurs through estrogen action

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use Biotécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiry date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product or any technical doubt on handling this product or this procedure,

contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotecnicald.com.br

AUTOMATI ON

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.biotecnicald.com.br

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