## C3 Plus

C3 Plus | C3 Plus  
Ref. 20.022.00

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### ENGLISH

#### INTENDED USE

Kit for the quantitative determination of Complement 3 in serum and plasma. In vitro diagnostic use.

#### STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

#### WORKING PRINCIPLE

Immunoturbidimetric method in which anti-human C3 complement antibodies form insoluble complexes when mixed with samples containing Complement C3. The turbidity intensity is proportional to the Complement C3 concentration of the sample.

#### SAMPLES – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

**Sample Type:** Serum and plasma (EDTA or heparin).

**Collection, handling and preparation:** Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

**Preservation:** Complement C3 in serum and plasma is stable for 8 days at a temperature of 4 to 8 °C or -20 °C. Do not use lipemic sample or with hemolysis.

#### PRODUCT DESCRIPTION

<b>R 1</b>	Anti-complement goat antibodies C3; Tris buffer pH 8.2; preservative.	
------------	---	--

#### STABILITY IN USE

- After opening, the product (R1) in use is stable until the expiration date on the label, following the recommended storage conditions (2 to 8 °C).
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

#### TREATMENT AND HANDLING OF THE PRODUCT

##### A) REAGENT PREPARATION

R1: Ready-to-use reagent.

##### B) OPERATING INTERVAL

The operating range of the product is 13,28 mg/dL at 450,0 mg/dL.

For higher values, dilute the sample with NaCl 150 mM (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

### QUALITY CONTROL

The use of controls should be routine practice in the laboratory. For Calibration and Internal Quality Control Laboratory it is recommended to use the calibrator and control below:

Multicalibrador Plus	<b>REF</b>	21.011.00
Medium Multiparameter Control		21.003.00

### TEST PROCEDURE, CALCULATION AND INTERPRETATION

#### A) TEST PROCEDURE

REACTION PARAMETERS	
<b>Type of reaction</b>	Final Point / Growing
<b>Wave length</b>	340 nm
<b>Temperature</b>	37 °C
<b>Sample volume/CAL*</b>	2,5 µL
<b>R1 volume*</b>	250 µL
<b>Incubation time</b>	2 minutes
<b>Calibration</b>	Calibration curve**

\* The sample volumes / CAL and R1 can be modified as long as the established ratio is maintained.

#### \*\*CALIBRATION CURVE:

Dilute the PLUS Multicalibrador in NaCl 9g/L as follows:

Dilution	1	2	3	4	5	6
<b>Calibrador (µL)</b>	-	10	25	50	75	100
<b>NaCl 9g/L (µL)</b>	100	90	75	50	25	-
<b>Fator</b>	0	0,1	0,25	0,5	0,75	1

Multiply the concentration of Complement C3 indicated in the instruction of use of the Multicalibrador PLUS by the corresponding factor to obtain the concentration of Complement C3 of each dilution.

#### B) CALCULATIONS

The equipment automatically calculates the analyte concentration in each sample.

#### C) INTERPRETATION

C3 is the complement component of higher concentration in plasma and participates in the classical and alternative pathways of activation of this system. It is a protein synthesized by the liver, although bacterial endotoxins induce its synthesis through monocytes and fibroblasts. Increases and decreases in concentration are of great clinical significance. Increased levels of C3 concentration are related to acute phase responses (trauma, inflammation), biliary obstruction, and focal glomerulosclerosis. Decreases in concentration are related to deficiencies of genetic origin (risk for bacterial infection, especially by encapsulated bacteria), or acquired (vascular disease and severe infections).

### INTERFERING AND LIMITATIONS

**Hemolysis, Jaundice and Lipemia:** Hemoglobin > 300 mg / dL / Bilirubin > 20 mg / dL. Lipemia interferes with the assay.

**Medications:** see recommended literature reference (Young, 2000).

### PERFORMANCE CHARACTERISTICS

**Sensitivity:** Limit of detection: 1,845 mg/dL / Limit of quantification: 13,288 mg/dL.

**Prozone High Dose Effect:** There is no prozone effect up to 680 mg / dL Complement C3.

**Analytical Specificity:** The product specifically determines Complement C3 in the presence of other interfering substances in the sample, up to the concentrations reported above.

**Accuracy:** The method was compared with a similar method by the determination of 40 samples in duplicate. The regression equation  $y = 0,988x + 4,002$  and correlation coefficient  $r = 0,9906$  was obtained. Using this equation the estimated systematic error of 2,301% for a level of 160 mg / dL and 1,134% for a level of 300 mg / dL.

**Accuracy:** It was determined using samples at 03 decimal levels, two runs in duplicate per day for 20 days, being obtained.

Samples (mg/dL)	Repetitions	In-run Precision		Total Precision	
		SD (mg/dL)	%CV	SD (mg/dL)	%CV
78,994	80	0,874	1,1	0,907	1,1
169,585	80	2,247	1,3	3,099	1,8
388,450	80	3,738	1,0	5,041	1,3

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

## RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the reagents out of the specified storage conditions when they are not in use.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator, and reagent.
- The water level of the water bath should be higher than that of the test tubes containing the reactions.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

## REFERENCE RANGES

Adults (20 - 60 years old)	90 - 180 mg/dL
----------------------------	----------------

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

## MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer for reading at 340 nm.
- Glass and / or automatic pipettes.
- Test tubes
- Multicalibrator Plus Biotécnica – CAT BT 21.011.00

## ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

## QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use Biotécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiry date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: [sac@biotechnicaltda.com.br](mailto:sac@biotechnicaltda.com.br)

## AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at [www.biotechnicaltda.ind.br](http://www.biotechnicaltda.ind.br)

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- LAMBRIS, J.D., *Immunology today*. v. 9, p. 387-393, 1988.
- PRICE, C. P.; et al. *Ann Clin Biochem*. v. 20, p. 1-14, 1983.
- DATI, F. et al; *Eur J Clin Chem Clin Biochem*. v. 34, p. 517-520, 1996.
- TIEZ, Textbook of Clinical Chemistry. 3 ed. Burtis CA, Ashwood ER. **WB Saunders** Co, 1999.
- FRIEDMAN AND YOUNG, Effects of the disease on clinical laboratory tests, 3 ed., **AACC Press**, 1997.
- YOUNG, D.S. *Effects of drugs on clinical laboratory tests - vol. 2, 5 ed.* Washington DC: **AACC Press**, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem*. v.27 p.493-501, 1981.

## TABELA DE SÍMBOLOS INTERNACIONAIS / TABLE OF INTERNATIONAL SYMBOLS / TABELA DE SÍMBOLOS INTERNACIONALES

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