**Fator Reumatóide de Turbílatex**  
RF Turbílatex / Fator Reumatóide de Turbílatex  
Ref. 20.007.00

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**ENGLISH**

**INTENDED USE**

Kit intended to determine of the Rheumatoid Factors (RF) in serum.

**PRINCIPLE OF THE METHOD**

Latex particles covered with Immunoglobulin G are agglutinated by the rheumatoid factors (RF) present in the sample. The agglutination produces increment of the absorbance and is proportional to the RF concentration in the sample and may be determined by comparison with a calibrator of known concentration at 650nm.

**REAGENT COMPOSITION**

<b>R</b> <b>DIL</b>	TRS buffer 20 mmol/L pH 8,2, preservatives	
<b>R</b> <b>LAX</b>	Latex particles covered with Human Immunoglobulin G pH 7.4, preservatives	
<b>CAL</b>	Human serum containing rheumatoid factors; preservatives	

All the components of human origin show negative results to HBs antigen, and anti-HV. However, it may be treated as potentially infectious material.

**STORAGE AND STABILITY**

- After first opening the reagents are stable until the expiration date printed on the label if stored 2-8°C tightly closed
- Do not use reagents over the expiration date
- Store at 2 to 8°C
- Reagents must remain outside the specified temperature only the time required for testing
- Do not freeze and protect from light.

**ADDITIONAL EQUIPMENT**

- Spectrophotometer or photometer with a thermostatic cuvet able to read at 650 (600-650) nm
- Thermostatic water bath at 37°C
- Pipettes and micropipettes.
- Clock or chronometer.
- Assay tubes.

**WARNINGS AND PRECAUTIONS**

- Perform the sample collection according to Good Laboratory Practices
- The kit is intended only for in vitro diagnostic use
- The samples must be considered as potentially infective.
- Wear the individual protective equipment according to Good Laboratory Practices
- Discard the reaction surplus equipment according to Good Laboratory Practices in a proper place for potentially infective material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product. The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, also available on the site [www.biotecnica.ind.br](http://www.biotecnica.ind.br) or calling for 55-35-3214-4646
- Do not combine reagents from different batches
- Do not exchange the caps from different reagents in order to avoid cross contamination
- Do not use the reagent if it displays any significant disagreement with the one specified in MSDS of the product
- Use glass pipettes and disposable tips for each sample, control, standard/calibrator and reagent.
- The level of water in the thermostatic bath must be over the level of mix reaction with the assay tubes.
- The laboratories must establish the chemical requirements, microbiological and particulate matter for water prior to and define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The appropriate cleaning and drying of the lab material are fundamental factors for the stability of the reagents to reach correct results.

**SAMPLE - PREPARATION AND STABILITY**

Serum  
Stable 7 days between 2 to 8°C or 3 months at -20°C  
Discard lipemic and/or hemolyzed samples.  
Samples containing fibrin must be centrifuged for 10 minutes at 3.000 rpm. Separate the supernatant and carry out the test.

**TECHNICAL PREPARATION**

**A) REAGENT PREPARATION**

The reagents (RD L e RLAX) are ready to use.  
Do not freeze the reagent. It must affect its stability and turn it unavailable for use.

**RF Calibrator**

The calibrator is lyophilized

1. Let the calibrator at room temperature (10 - 30°C) for 10 minutes.
2. Reconstitute with 2.0 mL of purified water.
3. Dissolve with gentle agitation. Let rest for 10 minutes, avoid frothing
4. Mix gently by inversion until the complete dissolution. The calibrator is stable for 30 days at 2 to 8°C and 3 months if stored at -20°C

**B) PROCEDURE**

**B1 - CURVE CALIBRATION**

Proceed with calibrator dilution of rheumatoid factor with saline solution (NaCl 0.9%), as follows:

Dilution	1	2	3	4	5	6
Calibrator (μL)	-	25	50	100	200	400