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ENGLISH
INTENDED USE
Kit intended to determine concentration of Ferritin in serum. Diagnostic use only.

- STORAGE AND HANDLING**
- Store at 2 to 8°C
 - Protect from light.
 - Stable until the kit expiration date that is printed on the label.
 - Do not use reagents whose date has expired.

WORKING PRINCIPLE
Latex particles covered with antibodies anti-human Ferritin are agglutinated by the Ferritin present in the sample. The agglutination produces increment of the absorbance and is proportional to the Ferritin concentration in the sample and may be determined by comparison with a calibrator of known concentration.

SAMPLE – PREPARATION AND STABILITY
Sample Type: Serum
Collection, handling and preparation: Make the collection of the sample in accordance with Good Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material. Discard lipemic and/or hemolyzed samples. Samples containing fibrin must be centrifuged before assay.
Preservation: stable in serum for 7 days from 2 to 8°C and 3 months from -20°C

PRODUCT DESCRIPTION	
 RDIL	TRS buffer ≥ 20 mmol/L pH 8,2; preservative
	

 RLAX	Latex particles suspended covered with antibodies anti-human Ferritin pH 8,2; preservative	
 CAL	Human serum containing Ferritin; preservative	

The determination of Urea is traceable to reference material 2nd International Ferritin Standard (85/578) from WHO (1992).

- STABILITY IN USE**
- The stability of the product (RDIL, RLAX e CAL) in use is 24 months as long as followed by the recommended storage conditions (2 to 8°C).
 - Reagents must remain outside the specified temperature only the time required for testing

TECHNICAL PREPARATION
A) REAGENT PREPARATION

Diluent and Latex (RDIL / RLAX)
Reagents ready for use (RDIL and RLAX). The freezing of the reagent affects the functionality of it, unavailable for use. The presence of small dark particles in the latex results from the manufacturing process and does not interfere in the product performance.

Ferritin Calibrator
The calibrator is lyophilized. Remove from the refrigerator and leave at 10-30°C for 10 minutes. Reconstitute with 3,0 mL of purified water. Dissolve with gentle shaking and leave it to stand for 10 minutes to prevent foam. Homogenize gently by inversion until completely dissolved. After reconstitution, the calibrator is stable for 30 days if stored at 2 to 8°C and 3 months if kept at -20°C

B) OPERATING RANGE
The product operating range is from 8,6 µg/L to 600 µg/L. For higher values, dilute the sample with NaCl 150 mM (0.9%), repeat the assay and multiply the result by the dilution factor.

QUALITY CONTROL
Use of controls should be a routine practice in the lab. It is suggested to use a control within the reference range or decision level and another control with value in another range of clinical significance. For the internal quality control of the lab it is indicated the use of the calibrator serum and control sera below

Multi-level Multi-parameter Control		21.003.00
-------------------------------------	---	-----------

TEST PROCEDURE, CALCULATION AND INTERPRETATION
A) TEST PROCEDURE

Calibration curve
Proceed dilution of Ferritin calibrator with saline solution (NaCl 0.9%), as follows:

Dilution	1	2	3	4	5	6
Calibrador (µL)	-	25	50	100	200	400
NaCl 9g/L (µL)	400	375	350	300	200	-
Factor	0	0,0625	0,125	0,25	0,5	1,0

For the concentration of each dilution of ferritin, dividing the concentration of the calibrator by the corresponding factor indicated in the table

- Test**
- Pre-heat the reagents to be used in the assay for one minute at 37°C
 - Zero the equipment with purified water.
 - Read the test tube.

RDIL	800 µL
RLAX	200 µL
CAL or Sample	90 µL

- Mix and immediately insert in the cuvette holder thermostatically at 37°C
- Measure the initial absorbance (A1) and at 5 minutes (A2) at 540 nm (530-550).

B) CALCULATIONS
Calculate the difference absorbance A2 - A1 (ΔA) of each calibrator dilution. Make a calibration curve by plotting on the X axis the concentration of each dilution and on the Y, the value of the respective ΔA. Calculate the difference in absorbance A2 - A1 (ΔA) of each sample. Determine the concentration of samples by interpolating the ΔA obtained from the calibration curve.

Q) INTERPRETATION
Serum ferritin concentration usually reflects body iron stores and is considered one of the most reliable indicators of iron status of patients. Whereas low serum concentrations of ferritin are always indicative of an iron deficiency, elevated concentrations can occur for a variety of reasons. Thus, although elevated concentration often indicate an excessive iron intake, they are also caused by liver disease, chronic inflammation and malignancies. Pregnant women, blood donors, hemodialysis patients, acid ascents and children are groups particularly at risk.

INTERFERING AND LIMITATIONS
Anti-coagulants: They cause measurement errors.