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ENGLISH

INTENDED USE
Kit intended to the quantitative determination of the Hemoglobin A1c (HbA1c) in human blood. Diagnostic use only.

STORAGE AND HANDLING
Store at 2 to 8 °C.
Protect from light.
Do not freeze.

Stable until the kit expiration date that is printed on the label.
Do not use reagents whose date has expired.

WORKING PRINCIPLE

This method is based in the interaction of antigen-antibody to the direct determination of HbA1c in total blood. The total hemoglobin has the same nonspecific absorption rate of the latex particles present in R1. When the mouse monoclonal antibody anti-human HbA1c is added (R2) the complex latex+HbA1c+mouse monoclonal antibody anti-human HbA1c is formed. The agglutination occurs when the goat polydonal antibody interacts with the mouse monoclonal antibody. Agglutination leads to turbidity which can be measured in absorbance spectrophotometer.

SAMPLE - PREPARATION AND STABILITY

Sample Type: Total Blood with EDTA

Collection, handling and preparation: Make the collection of the sample in accordance with Good Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: Stable for 7 days at 2 to 8 °C

PRODUCT DESCRIPTION

R 1	Latex, glycol buffer.	
R 2	Glycol buffer, Mouse monoclonal antibody anti-human HbA1c, goat polydonal antibody anti-mouse IgG stabilizers.	
R HML	Purified Water, Stabilizer agents.	
CONTROL	Human total blood, glucose, stabilizer agents.	

STABILITY USE

- R1, R2 and CONTROL in use is 30 days as long as followed by the recommended storage conditions (2 to 8 °C).
- RHM in use is 24 months as long as followed by the recommended storage conditions (2 to 8 °C).

- Reagents must remain outside the specified temperature only the time required for testing.

TECHNICAL PROCEDURE

A) REAGENT PREPARATION
R1 and R2
R1 and R2 are ready to use. Homogenize gently.

CONTROL

The control is lyophilized. Put at room temperature (10 to 30 °C) for 10 minutes. Reconstitute the lyophilized with 0,5 mL of purified water. Let it rest for 30 minutes. Homogenize it gently by inversion until the full dissolution. After reconstitution the control is stable for 30 days if stored at 2 to 8 °C

B) OPERATING RANGE

The product operating range is from 2% to 16%

QUALITY CONTROL

Use of controls should be a routine practice in the lab. It is suggested to use a control within the reference range or decision level and another control with value in another range of clinical significance. For the internal quality control of the lab it is indicated the use of the calibrator serum and control sera below.

HbA1c Calibrator Set	21.009.00
HbA1c Plus CONTROL	20.025.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

HEMOLYSIS

- The procedure of hemolysis must be applied on samples, control and calibrators, as described below:
- 1 - Pipette 1,0 mL of RHM in assay tubes to each sample, control and calibrator.
 - 2 - Homogenize the Total Blood, pipette 20 µL in each tube and shake.
 - 3 - Let the tubes rest for 5 minutes to complete the hemolysis fully. The hemolyzed samples are stable for 7 days at 2 to 8 °C

CALIBRATION CURVE

Use the HbA1c Calibrator on Set. CAT BT 21.009.00

Point 0	NaCl 0,9%
Point 1	Level 1
Point 2	Level 2
Point 3	Level 3
Point 4	Level 4

Note: It's necessary to proceed with calibrators hemolysis as described above.

ASSAY

- 1- Heat the Reagent 1 for 3 minutes at 37 °C
- 2- Pipette in the assay tube:

	Blank	Calibrators	Sample
R1	720 µL	720 µL	720 µL
Calibrator	---	20 µL	---
Sample	---	---	20 µL

- 3- Incubate 5 minutes at 37 °C

Reagent 2	240 µL	240 µL	240 µL

4- Homogenize and incubate 5 minutes at 37 °C
5- Measure the absorbance of the calibrators (in increasing order of concentration) and of the samples against the blank at 630 nm (600-660 nm).

B) CALCULATION
Obtain the calibration curve by plotting on X axis the concentration of each calibrator and on Y axis the value of the respective absorbance. Determine the concentration of the samples by interpretation of the absorbance sample values in the calibration curve.

Q) INTERPRETATION

Glycohemoglobin is the most important marker to monitor diabetic patients. During the circulatory life of the red cells the HbA1c is continually formed by the addition of glucose in the non-terminal portion of the beta-subunit of hemoglobin. This process reflects an average of the exposition of the hemoglobin to the glucose for a long period. Glycohemoglobin formation is dependent of glucose levels and half-life of hematis (about 120 days). Thus it reflects the glucose levels of 6 to 8 weeks before the dosage. This dosage must be used exclusively to monitor diabetic patients, and can not be used in to confirm the diabetes mellitus diagnosis. Factor that influence the half-life of hematis as well as hemoglobinopatias are possible interferences to the

dosage This method uses the antigen – antibody interaction, determining exclusively the HbA1c fraction.

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Drugs: see recommended reference (Young, 2000).

Inconsistent results may occur in patients in the following conditions:

- Use of opiates, diuretics, alcoholism, high doses of acetylsalicylic acid
- High levels of fat in hemoglobin (HbF) may lead to decreased results of HbA1c
- The variant hemoglobin: Hemoglobin S (HbS) and Hemoglobin A2 (HbA2) and labile intermed品 are from glycosylated hemoglobin (Schiff base) are not detected in this immunoassay and do not interfere in the determination.

PERFORMANCE CHARACTERISTICS

Sensitivity: 1,8%
Analytical Specificity: The product determines HbA1c specifically in the presence of other interfering substances in the sample up to the concentrations stated above.

Accuracy: The method was compared with a similar method for the determination of 45 samples in duplicate. The regression equation obtained was y = 1,0695x - 0,311 and the correlation coefficient r = 0,9677.

Intra-Assay: The realization of 20 determinations of the same sample at the same day showed a Coefficient of Variation of 3,83%.

Inter-Assay: The realization of 10 determinations of the same sample at different days showed a Coefficient of Variation of 1,82%.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with Good Laboratory Practice
- Do not mix reagents from different lots.
- Do not exchange the caps from different reagents in order to avoid cross contamination
- Do not use the reagent if it displays any signs of damage with the one specified in the product MSDS.
- Avoid leaving reagents out of the specified storage conditions, when they are not in use.
- Use glass pipettes and disposable tips for each sample control, standard / calibrator and reagent.
- The level of the water bath must be greater than the test tubes containing the reaction.
- The level of water in the thermostatic bath must be over the level of mix reaction with in the assay tubes.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

REFERENCE RANGES

Non-diabetic individuals	4,0% to 6,0%
Diabetic individuals in Glucemic Control	
Adults	Less than 7,0%
0 up to 6 years old	Between 7,5 and 8,5%
6 up to 12 years old	Less than 8,0%
13 up to 19 years old	Less than 7,5%
Bilirubin	Less than 8,0%

These values are intended for orientation only. It is recommended that each lab establishes its own reference range.

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer with a thermostatic cuvette holder to read at 630 nm (600 – 660).
- Thermostatic water bath at 37 °C
- Glass pipettes and/or automatic
- Clock or chronometer.
- Assay tubes

WARNINGS AND PRECAUTIONS

- Discard the reaction surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.bioteccaltda.com or calling for 55-35-3214-4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use, biotechnical reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiry date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@bioteccaltda.com

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.bioteccaltda.com

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