IgG PLUS

IgG PLUS | IgG PLUS
Ref. 20.019.00

FINALIDADE

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ENGLISH

INTENDED USE

Kit intended to determination of Immunoglobulin G in serum and Plasma.

PRINCIPLE OF THE METHOD

The determination is based on a turbidimetric method in which anti-immunoglobulin G antibodies form with IgG an insoluble complex, providing a turbidity, whose intensity is proportional to the IgG concentration of the sample determined spectrophotometrically at 540 nm.

REAGENT COMPOSITION

R 1	Goat antibodies anti-Immunoglobulin G; Tris buffer pH 8.2; preservatives.	
------------	---	---

STORAGE AND STABILITY

- After first opening the reagents are stable until the expiration date printed on the label if stored 2-8 °C tightly closed.
- Do not use reagents over the expiration date.
- Store at 2 to 8 °C.
- Do not freeze and protect from light.
- Reagent must remain outside the specified temperature only the necessary time for the tests

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer able to read at 540 nm. (520-560).
- Thermostatic water bath at 37 °C.
- Pipettes and micropipettes.
- Clock or chronometer.
- Assay tubes.
- Multi Calibrator PLUS Biotécnica - CAT BT 21.011.00.

WARNINGS AND PRECAUTIONS

- Perform sample collection in accordance with Good Laboratory Practice
- The kit is intended only for in vitro diagnostic use.
- The samples must be considered as potentially infective.
- Wear the individual protection equipment accordingly to Good Laboratory Practice.
- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product. The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, also is available on the site www.biotecnica.ind.br or calling for 55-35-3214-4646.
- Do not combine reagents from different batches
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any sign in disagreement with the one specified in MSDS of the product
- Avoid leaving reagents outside the specified storage conditions, when they are not in use.
- Use glass pipettes and disposable tips for each sample, control, standard/calibrator and reagent.
- The level of water in the thermostatic bath must be over the level of mix reaction within the assay tubes.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

SAMPLE – PREPARATION AND STABILITY

Fresh serum or plasma with heparin or EDTA
Stable 7 days at 2 to 8 °C or 3 months at -20 °C.

Discard lipemic and/or hemolyzed samples.

Samples containing fibrin must be centrifuged for 10 minutes at 3,000 rpm. Separating the supernatant and conduct the assay.

TECHNICAL PROCEDURE

A) REAGENT PREPARATION

Ready-to-Use reagent

Note: The presence of particles, turbidity or absorbance of higher reagent over 0.300 at 540 nm, zeroed with purified water indicates deterioration.

B) CALIBRATION CURVE

Dilute Multicalibrador PLUS in NaCl 9 g / L as follows:

Dilution	1	2	3	4	5	6
Calibrador (µL)	-	10	25	50	75	100
NaCl 9g/L (µL)	100	90	75	50	25	-
Fator	0	0,1	0,25	0,5	0,75	1

Multiply the concentration of IgG Multicalibrador PLUS by the corresponding factor to obtain the concentration of each dilution of IgG.

C) PROCEDURE

- Heat R1 and the spectrophotometer (cuvet holder) at 37 °C.
- Adjust the blank (zero) of the equipment with deionized water at 540 nm.
- Pipette in the assay tube:

	Volume
Sample/ Calibrator	7 µL
R1	1,0 mL

- Mix and insert into the cuvette holder thermostatically.
- After 2 minutes, measuring absorbance (A) of the sample / calibrator.

6. The concentration of IgG in the sample is calculated by interpolation of their absorbance on the calibration curve.

CALIBRATION/CALCULATION PROCEDURES

- a) Make a calibration curve by plotting the X-axis the concentration of each calibrator dilution and the Y-axis value of the corresponding absorbance (A).
 b) Determining the concentration of samples by interpolating the absorbance readings (A) obtained for the samples from the calibration curve.

SENSITIVITY AND LINEARITY

Analytical Sensitivity: 3,25 mg/dL
 Linearity Limit: 3000 mg/dL

For higher values, dilute again the sample with NaCl solution (0.9%), repeat the assay and multiply the obtained result by the dilution factor.
 Prozone: The reaction has no prozone effect up to 7000 mg/dL.

TECHNICAL LIMITATIONS

Interferences:

- Anticoagulantes: They suffer interference from citrate, oxalate and fluoride.
- Bilirrubin > 40 mg/dL; Hemoglobin > 8g/L ; Triglycerides > 2000 mg/dL

QUALITY CONTROL

The use of controls to evaluate the imprecision of the analysis must be a routine practice in the lab. It is suggested to use a control within the reference range or decision level and a control with value in another range of clinical significance.

For Internal Calibration and Control it is recommended the use of the calibrator serum and control sera below:

Multi calibrator PLUS	REF	21.011.00
Multiparameter Control		21.003.00

REFERENCE RANGES

IgG Reference Values (mg/dL)	
Newborns	299 up to 852
Adults	700 up to 1600

These values are intended for guidance, recommended that each laboratory establish its own reference range.

PERFORMANCE CHARACTERISTICS

Intra-Assay

The realization of 20 determinations of the same sample at the same day showed a Coefficient of Variation of 2.50%.

Inter-Assay

The realization of 10 determinations of the same sample at different days showed a Coefficient of Variation of 3.46%.

Analytical Specificity

A comparison with a reference method showed a correlation coefficient (r) of 0.9636 obtained from ambulatory samples. The result equation of the linear regression is $Y = 1.005 X - 2,94$.

CLINICAL SIGNIFICANCE

The quantification of immunoglobulins in serum is important for the diagnosis of primary or secondary immunodeficiency, for therapy monitoring with immunoglobulin and to monitor the clinical course of multiple myeloma.

IgG is the most important immunoglobulin and represents about 75% of total circulating immunoglobulins.

The decrease in IgG levels occurs mainly acquired and congenital immunodeficiencies.

The elevation of polyclonal IgG (normal response to infection) occurs in hepatitis, cirrhosis, and autoimmune diseases. Increases of monoclonal IgG (paraprotein) are found in multiple myeloma and lymphocytic leukemia.

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotecnicaltda.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.biotecnicaltda.ind.br

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