# IgM PLUS

## IgM PLUS | IgM PLUS

Ref. 20.020.00

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#### ENGLISH

#### INTENDED USE

Kit for the determination of Immunoglobulin M (IgM) in serum and plasma. In vitro diagnostic use.

#### STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

#### WORKING PRINCIPLE

Immunoturbidimetric method wherein human anti-immunoglobulin M antibodies form insoluble complexes when mixed with samples containing IgM. The turbidity intensity is proportional to the IgM concentration of the sample.

#### SAMPLES – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

**Sample Type:** Serum and plasma (heparin or EDTA).

**Collection, handling and preparation:** Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

**Preservation:** IgM in serum and plasma is stable for 4 months at 4 to 8 ° C and 6 months at -20 ° C. Do not use lipemic sample

#### PRODUCT DESCRIPTION

<b>R 1</b>	Goat anti-immunoglobulin M antibodies; Tris buffer 20 mmol / L pH 8.2; preservative.	
------------	--	---

#### STABILITY IN USE

- After opening, the product (R1) in use is stable until the expiration date on the label, following the recommended storage conditions (2 to 8 ° C).
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

#### TREATMENT AND HANDLING OF THE PRODUCT

##### A) REAGENT PREPARATION

R1: Ready-to-use reagent.

##### B) OPERATING INTERVAL

The operating range of the product is from 14.25 mg / dL to 399.21 mg / dL. For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

#### QUALITY CONTROL

The use of controls should be routine practice in the laboratory. For Calibration and Internal Quality Control Laboratory it is recommended to use the calibrator and control below:

Multicalibrador Plus	<b>REF</b>	21.011.00
Medium multicalibrador Plus		21.003.00

#### TEST PROCEDURE, CALCULATION AND INTERPRETATION

##### A) TEST PROCEDURE

REACTION PARAMETERS	
Reaction Type	Final/Growing Point
Wave length	340 nm
Temperature	37 °C
Sample Volume/CAL*	2 µL
R1* Volume	200 µL
Incubation Time	2 minutes
Calibration	Calibration Curve**

\* Sample volumes / CAL and R1 may be modified since established ratio is maintained.

##### \*\* CALIBRATION CURVE:

Dilute the PLUS Multicalibrador in 9g / L NaCl as follows:

Dilution	1	2	3	4	5	6
Calibrador (µL)	-	10	25	50	75	100
NaCl 9g/L (µL)	100	90	75	50	25	-
Factor	0	0,1	0,25	0,5	0,75	1

Multiply the IgM concentration indicated in the Multicalibrador PLUS use instruction by the corresponding factor to obtain the IgM concentration of each dilution. The IgM concentration of the calibrator is traceable to Standard Reference Materials ERM DA470k / IFCC.

##### B) CALCULATIONS

The equipment automatically calculates the analyte concentration in each sample.

##### C) INTERPRETATION

Quantification of serum immunoglobulins is important for the diagnosis of primary or secondary immunodeficiencies, for the monitoring of immunoglobulin therapy, and for the follow-up of the clinical course of multiple myeloma. IgM is the only immuno-globulin synthesized by newborns. It accounts for about 5% to 10% of total circulating immunoglobulins. It is a pentameric molecule, and its large size prevents its passage into extravascular spaces. Its decrease occurs in acquired or congenital immunodeficiencies. Elevation of polyclonal IgM (normal response to infections) occurs primarily in primary viral infections and infections of the bloodstream, such as malaria. It may also be increased in primary biliary cirrhosis and chronic active hepatitis. Increases in monoclonal IgM (paraprotein) are found in plasma cell proliferative disorders, such as Waldenström macroglobulinemia.

#### INTERFERING AND LIMITATIONS

**Hemolysis, Jaundice and Lipemia:** Hemoglobin > 200 mg / dL, Bilirubin > 20 mg / dL and lipemia interfere with the test.

**Medications:** consult recommended literature reference (Young, 2000).

#### PERFORMANCE CHARACTERISTICS

**Sensitivity:** Detection limit: 2,424 mg / dL / Limit of quantification: 7,685 mg / dL. **Prozone High Dose Effect:** Prozone effect does not occur up to 747,24 mg / dL IgM.

**Analytical Specificity:** The product specifically determines IgM in the presence of other interfering substances in the sample up to the concentrations reported above.

**Accuracy:** The method was compared with a similar method by the determination of 23 samples in duplicate. The regression equation  $y = 0,996x - 0,50$  and correlation coefficient  $r = 0,9755$  was obtained. Using this equation the estimated total systematic error of -0.73% for a level of 150 mg / dL and - 0.52% for a level of 400 mg / dL.

**Accuracy:** It was determined using samples at 03 decimal levels, two runs in duplicate per day for 20 days, being obtained.

Samples (mg/dL)	Repetitions	In-run Precision		Total Precision	
		SD (mg/dL)	%CV	SD (mg/L)	%CV
86,086	80	1,097	1,3	1,758	2,0
175,392	80	2,492	1,4	3,715	2,1
316,677	80	0,834	0,3	3,398	1,1

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

#### RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the reagents out of the specified storage conditions when they are not in use.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

#### REFERENCE RANGES

Newborns ( 4 days )	5 to 30 mg/dL
Adults (20-60 years)	40 to 230 mg/dL
Adults ( > 60 years)	30 to 360 mg/dL

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

#### MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer for reading in 340.
- Glass and / or automatic pipettes.
- Test tubes.
- Biotécnica PLUS Multicalibrador - CAT BT 21.011.00

#### ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at [www.biotecnica.ind.br](http://www.biotecnica.ind.br) or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: [sac@biotecnicaltda.com.br](mailto:sac@biotecnicaltda.com.br)

#### AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at [www.biotecnicaltda.ind.br](http://www.biotecnicaltda.ind.br)

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

- BURTIS, C.A.; ASHWOOD, E.R.; BRUNS, D.E. Tietz: Fundamentos de Química Clínica. 6. ed. Rio de Janeiro: Saunders Elsevier, 2008. 959 p.
- DATI, F. *et al*; Consensus of a Group of Professional Societies and Diagnostic Companies on Guidelines for Interim Reference Ranges for 14 Proteins in Serum Based on the Standardization against the IFCC/BCR/CAP Reference Material (CRM 470). *Eur J Clin Chem Clin Biochem.* 34, 517-520 (1996).
- PONTET F.; AUGET, J.L.; ROUSSELET, F. Six Approaches to Routine Quantification of Human Monoclonal Immunoglobulin M Compared. *Clin. Chem.*, 30/5. 774-777 (1984).
- SKOUG, J.W.; PARDUE, H.L. Kinetic Turbidimetric Method for the Immunochemical Quantification of Immunoglobulins, Including Samples with Excess Antigen. *Clin. Chem.* 34/2, 309-315 (1988).
- WESTGARD, J. O. *et al*. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACC Press, 2000.

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