Afa Amilase Responsável Técnico
Dr. Glison Séri Rizzo
CRF MG – 5310
MS 80027310202
Alpha Amylase / Alfa Amilase
Ref. 11.001.00

R NALIDADE

Kit destinado a determinar a atividade enzimática da alfa-amilase no soro, plasma, urina, líquido do ascítico e Pleural.

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ENGLISH

INTENDED USE

Kit intended to determine of Alpha amylase activity in serum plasma, urine, pleural and ascitic fluids.

PRINCIPLE OF THE METHOD

This kit uses the substrate G-Glú-CNIP (2-chloro-4-nitrophenyl-galactopyranosil maltodose) which is hydrolyzed by the Alpha-amylase from sample releasing the Cl or 4-nitrophenol (CNP). The CNP absorbs light at 405 nm and the absorbance increase is measured as the enzymatic activity.

REAGENT COMPOSITION

R 1	MES buffer 50 mmol/L pH 6,00, Sodium Chloride 70 mmol/L; Calcium Chloride 6 mmol/L; Gal-Glú-CNIP 2,22 mmol/L; Sodium Azide 0,095%	X
------------	---	----------

STORAGE AND STABILITY

- The reagents are stable up to the date stated on the label.
- Do not use reagents over the expiration date.
- Store at 2 to 8 °C
- Do not freeze and protect from light.

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer with a thermostatic cuvette able to read at 405 nm

- Water bath, thermostatically controlled at 37 °C.
- Pipettes and micropipettes.
- Clock or chronometer.
- The reagents must remain outside of the specified temperature only the necessary time to carry out the tests.
- Assay tubes.

WARNINGS AND PRECAUTIONS

- The kit is intended for in vitro diagnostic use only.
- All speed mes and controls should be handled as potentially infective.
- Wear the individual protective equipment according to Good Laboratory Practice.
- Discard the reaction surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for Disposing Security and first Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotecnica.ind.br or calling for 55-35-3214-4646.
- Do not combine reagents from different batches.
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any sign of contamination.
- The reagents must be kept out of specified temperatures only the time necessary for the realization of the tests.
- Use glass pipettes and disposable tips for each sample control, standard/calibrator and reagent.
- The level of water in the thermostatic bath must be over the level of mix reaction with in the assay tubes.
- Precautions to observe in order to avoid the contamination of Ri: Do not pipette the reagent by mouth, do not wipe, do not use glassware contaminated with saliva or sweat. Do not talk near the uncapped flasks since it may contaminate reversibly the reagents. It is advisable to use latex gloves. After using the reagent cap immediately and store at 2 to 8 °C protected from the light.

SAMPLE - PREPARATION AND STABILITY

- Serum or plasma (Heparin): The alpha-amylase in serum and plasma (heparin) is stable for 7 days at 2 to 8 °C
- Urine: Work with samples collected within the period of 24 hours. It's not necessary to add preservatives.
- Centrifuge the sample for 10 minutes at 3000 rpm
- Execute the assay using the supernatant.

The alpha-amylase in Urine is stable for 7 days at 2 to 8 °C at 7,0 pH. Adjust the pH with HCl 0,1 N or NaOH 0,1 sd utions.

- Ascitic and Pleural Fluids
- Centrifuge the sample for 10 minutes at 3000 rpm
- Execute the assay using the supernatant.

TECHNICAL PREPARATION

A) REAGENT PREPARATION

The reagent is ready to use
If closed, the reagent is stable up to expiration date stated on the label, if stored at 2 to 8 °C. Once opened the reagent is stable 30 days, since it is designed immediately after use and stored at 2 to 8 °C. During handling the reagent is susceptible to chemical and microbiological contamination which may reduce the reagent stability.

Note: The absorbance of the reagent above 0,600 at 405 nm with the instrument adjusted to zero with still red or deionized water, indicates deterioration.

B) PROCEDURE

1. Heat the reagent for 3 minutes at 37 °C
2. Pipette in the assay tube:

R 1	Volumen
	1,0 mL
Sample	10 µL

3. Homogeneizar e inserir imediatamente no termostaticado cuvette
4. Após 1 min, note the initial absorbance (A) and read again after exactly 1, 2 and 3 minutes (A, A e A, respectively).

CALIBRATION/CALCULATION PROCEDURES

Using the measured absorbances, calculate the mean variation of absorbance per minute ($\Delta A/n\text{min}$):

$$(\Delta A/n\text{min}) = \frac{(A - A) + (A - A) + (A - A)}{3}$$

The activity of the Alpha-amylase is then calculated by multiplying the $\Delta A/n\text{min}$ by the factor below

$$\text{Alpha-amylase (U/L)} = \Delta A/n\text{min} \times 6196$$

Uridin Amylase per Hour
Uridin Amylase/Hour (U/H) = Amylase (U/L) x urinary volumen (mL) / 1000 x collection time (hours)

Note: To calculate the Factor for automatic analyzers use the absorptivity coefficient of Gal-Glú-CNIP = 16,3

SENSITIVITY AND LINEARITY

- Methodology of Sensitivity: Cálculo do método: 6,19 U/L / Experimental method: 3,167 U/L - Linearity: 1000 U/L. For higher values, use the sample with NaCl 150 mM (0,9%) sd ution, repeat the assay and multiply the obtained result by the dilution factor.

TECHNICAL LIMITATIONS

- Anticoagulants: Fluoride, oxalate, citrate and EDTA may interfere with the Alpha-amylase activity leading to false low results.
- Hemolyzed, Jaundiced and Lipemic Serum: Hemoglobina ≥ 180 mg/dL / Bilirrubina ≥ 20 mg/dL / Triglicérides ≥ 1500 mg/dL
- Medication Effects: False Elevated Results: Meprednolona, codeína, morfina, álcool, diuréticos, sais dilúos e tetracyclina

QUALITY CONTRL	
Any clinical laboratory must keep an internal quality control program which defines the aids, procedures and criteria for the tolerance limits, corrective actions and registration of the activities. Also it must be kept a defined system of verification of the analytical variability that occurs in any measuring system. The use of controls to evaluate the imprecision of the analysis must be a routine practice in the lab. It's suggested to use a control within the reference range and other control within the clinical significance range. The application of the Multiple Rules of Westgard for evaluation of the control state is recommended. The lab must participate of external quality control programs. For the internal quality control of the lab it's indicated the use of the calibrator serum control sera below	
Calibrator Serum- Autoal H	13.002.00
Normal Control Serum- Quanti norm	REF 13.003.00
Pathological Control Serum- Quanti At	13.004.00

REFERENCE RANGES	
	ACTIVITY AT 37 °C
Serum	< 86 U/L
Urine	< 470 U/L
Ascitic Fluid	20 a 104 U/L
Pericardial Fluid	20 a 104 U/L

These values are intended for orientation only. It's recommended that each lab establishes its own reference range.

Conversion for International System of Units (SI): $\mu\text{kat/L}$
 A-pha-amylase (U/L) \times Q.017 = A-pha-amylase ($\mu\text{kat/L}$)

PERFORMANCE CHARACTERISTICS

Intra-Assay
 The measurement of 20 determinations of the same sample at the same day showed a Coefficient of Variation of 3.81%

Inter-Assay
 The measurement of 10 determinations of the same sample at different days showed a Coefficient of Variation of 2,04%

Analytical Specificity
 A comparison with a reference method showed a correlation coefficient (r) of 0.995 obtained from analytical samples. The result equation of the linear regression is $Y = 1,0031 X - 3,3254$.

CLINICAL SIGNIFICANCE

The A-pha-amylase is an enzyme that hydrolyzes carbohydrate complexes. There are two A-pha-amylase isoenzymes being one from pancreatic and other from salivary gland origin. It's known that pancreatic and salivary A-pha-amylases are present in a 40:60 proportion in healthy individuals. The determination of serum and urinary A-pha-amylases are widely used in the diagnosis of pancreas diseases as well as investigations of pancreatic function. In the most patients with acute pancreas inflammation, the serum level of A-pha-amylase increases from 2 to 12 hours after the beginning of the inflammation. It reaches a peak in 24 hours and come back to normal values between 48 and 72 hours. The magnitude of the A-pha-amylase serum increases not related to the gravity of the pancreas inflammation but indicates, with high probability, the diagnosis of acute pancreas inflammation. However, the increase of serum A-pha-amylase is not necessarily due to acute pancreas inflammation. Lung and ovarian tumors may cause hyperamylase activity up to 50 times the reference value.

It's estimated that approximately 25% of the serum A-pha-amylase is eliminated in the urine. In kidney injury the A-pha-amylase activity is increased according to the extension of the kidney lesion. Normal A-pha-amylase serum levels may be found in patients with acute pancreas inflammation. In chronic pancreas inflammation the serum A-pha-amylase levels are frequently normal. The A-pha-amylase may bind to serum proteins (such as immunoglobulin) forming high molecular weight complexes known as macroamylases. This is characterized as high serum levels and normal urinary levels and may occur in healthy individuals.

Increased Levels	Decreased levels
Parotid inflammation	Pancreatic insufficiency
Acute pancreas inflammation	Advanced Cystic Fibrosis
Macroamylasemia	Serious Liver injury
Ectopic pregnancy	
Mesenteric occlusion	
Serious burns	
Salivary gland lesion	
Intra-abdominal diseases	
Alcohol intoxication	
Serious kidney insufficiency	
Obstruction of gall bladder ves	
Intestinal obstruction	
Pancreatic duct obstruction	
Pancreas cancer	
Diabetic ketoacidosis	
Lung and ovarian neoplasy	

Urinary A-pha-Amlyase
 A high quantity of serum A-pha-amylase is excreted in the urine. Then, serum increase causes the urine increase. However, the urinary A-pha-amylase may reach more elevated levels and is maintained for longer periods if compared to the serum A-pha-amylase. The macroamylasemia is characterized as high serum levels and normal urine levels. The relation between the A-pha-amylase/creatinine "clearance" helps in the differential diagnosis of macroamylasemia (very low "clearance"). The normal relations is 1% to 4%. This relation is increased in acute pancreas inflammation, decreased in macroamylasemia and normal or discretely decreased in salivary hyperamylasemia.

OBSERVATIONS

B IOTÉCNICA INDCOM LTDA - CNPJ: 025340690001/20. Avenida Vespertino R. 600 - Indústria Miguel de Luca CEP 37072-030 Varigües - MG BRAS L Tel/fax +55 35 3214 4646 www.bioteccalind.com.br

1. The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.

2. The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

PRESENTATIONS

QUALITY ASSURANCE/CUSTOMER TECHNICAL SERVICE
 Before being approved for use, Biotécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiry date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@bioteccalind.com.br.

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.bioteccalind.com.br.

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- WINDEEN E.S. et al. Development of a direct assay for amylase. *Qin Chem*, v.34, n.10, p.2005-2008, 1988.
- RAUSCHER E. et al. Optimal mixed cations for determining activity concentration of alpha-amylase in serum. *J Clin Lab Diagn* 4-nitrophenyl malt heptasoside as substrate. *Qin Chem*, v.31, n.1 p.14-19, 1985
- YOUNG D.S. *Effects of drugs on clinical laboratory tests - vol. 2*. 5ed. Washington DC: AACCPress, 2000.
- WESTGARD J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Qin Chem*, v.27, p.493-501, 1981.

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