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ENGLISH

INTENDED USE
Kit intended to determine iron binding capacity in serum.

PRINCIPLE OF THE METHOD
The iron (Fe) is transported in the blood by a specific beta-globulin called Transferrin. Each Transferrin molecule is able to transport two iron (Fe³⁺) ions and in physiological conditions the Transferrin is around 30% saturated. In order to determine the iron binding capacity of the Transferrin molecule, the test sample is incubated mixed to an iron standard. These free iron ions will saturate the available sites of the Transferrin and the iron excess (not bound) will be determined by a colored complex with ferrozine, indicating the iron binding capacity of the patient.

REAGENT COMPOSITION

R 1	TR S buffer 25 mmol/L pH 8,50; phenoxetanol Q.1%	X
R 2	Ferozine 30 mmol/L; hydroxylamine 300 mmol/L	X
STD	Hydroxylamine 400 mmol/L Iron sulphate equivalent to 500 µg/dL of iron.	X

The concentration of this standard was determined using the NST 937 International Standard.

STORAGE AND STABILITY

- The reagents are stable up to the date stated on the label.
- Do not use reagents over the expiration date.
- The reagents must remain outside the specified temperature only the necessary time for the tests.
- Store at 15 to 30 °C
- Do not freeze and protect from light.

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer able to read at 560 (540-580 nm).
- Glass Pipettes and/or automatic.
- Water bath, thermostatically controlled at 37 °C
- Clock or chronometer.
- Assay tubes.

WARNINGS AND PRECAUTIONS

- The kit is intended only for in vitro diagnostic use.
- The safety procedures for handling biochemical materials must be observed.
- The samples must be considered as potentially infective.
- Wear the individual protection equipment according to Good Laboratory Practice.
- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotecnica.ind.br or calling for 55-35-3214-4646.
- In case of accidental leak and eye or skin contact wash with plenty of water.
- In case of ingestion seek for medical advice.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product.
- Do not combine reagents from different lots.
- Use glass pipettes and disposable tips for each sample, control, standard/calibrator and reagent.
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any sign of contamination.
- **All the glass ware must be free of iron ions. In order to reach that it must to submerge the glass ware for 6 hours in a nitric acid 10-15% or chloric acid at 0,1 N. And to eliminate the activity by rinsing with several times with purified water. Dry the glass ware in stainless steel or plastic supports.**

SAMPLE - PREPARATION AND STABILITY

- Serum

The iron binding capacity is stable for 3 days from 2 to 8 °C.
Avoid hemolyzed samples.

TECHNICAL PROCEDURE

A) REAGENT PREPARATION
The reagents are ready to use.

B) PROCEDURE

1. Pipette in the assay tubes:

	Blank	Standard	Sample
Reagent 1	1,0 mL	----	1,5 mL
Deionized water	1,5 mL	2,0 mL	----
Sample	----	----	0,5 mL
Standard	----	0,5 mL	0,5 mL

2. Homogenize gently and incubate 10 minutes at 37 °C
3. Measure the absorbance (A1) against the Blank at 560 nm (540-580 nm) matching the zero with the Blank.
4. Add:

	Blank	Standard	Sample
Reagent 2	0,1 mL	0,1 mL	0,1 mL

5. Homogenize gently and incubate 10 minutes at 37 °C
6. Measure the absorbance (A2) against the Blank at 560 nm (540-570 nm) matching the zero with the Blank.

CALIBRATION/CALCULATION PROCEDURES

For the calculations below must be determined the Iron obtained by the kit:
IRON FERROQ NE BT 12.005.00

U BC = Latent Iron Binding Capacity
T BC = Total Iron Binding Capacity
TSI = Transferrin Saturation Index

U BC (µg/dL) = 500 - [(A2 - A1) / (A2 - A1) X 500] Standard Absorbance

T BC (µg/dL) = U BC + serum iron

TSI (%) = serum iron X 100 / T BC

Examp le
A1 sample = Q 042 Standard Absorbance = Q 0.480
A2 sample = Q 391 Serum iron = 110 µg/dL

U BC (µg/dL) = 500 - [(Q 391 - Q 042) / Q 480 X 500] = 136 µg/dL

T BC (µg/dL) = 110 + 136 = 246 µg/dL

TSI (%) = 110 x 100 = 44.7% / 246

SENSITIVITY AND LINEARITY

Method optical sensitivity (U BC)
Calculated: 1.075 µg/ D
Linearity: 450 µg/dL

For values higher, dilute the sample with Q9% NaD 150mM solution, repeat the assay and multiply the obtained result by the dilution factor.

TECHNICAL LIMITATIONS

- Anti coagulants
All the anti coagulants interfere in the reaction
- Jaundice and Lipemic Sera
Total Bilirubin > 10 mg/dL
Lipemic > 900 mg/dL

QUALITY CONTROL

Any clinical laboratory must keep an internal quality control program which defines the aids, procedures and criteria for the tolerance limits, corrective actions and registration of the activities. Asq it must be kept a defined system for verification of the analytical variability that occurs in any measuring system
The use of controls to evaluate the imprecision of the analysis must be a routine practice in the lab. It's suggested to use a control with in the reference range and other control within the clinical significance range. The application of the Multi-Phase Rules of Westgard for evaluation of the control status is recommended.

The lab must participate of external quality control programs.
For the internal quality control of the laboratory indicated the use of the calibrated serum and control sera below

Normal Control Serum- Quantiron m **REF** 13.003.00
Pathological Control Serum- QuantiAt **REF** 13.004.00

REFERENCE RANGES

U BC	140 to 280 µg/dL
T BC	250 to 410 µg/dL
TSI	20 to 50 %

These values are intended for orientation only. It's recommended that each lab establish its own reference range

Conversion for International System of Units (SI): mmol/L
T BC (µg/dL) X Q 179 = mmol/L

PERFORMANCE CHARACTERISTICS

Intra-Assay
The realization of 20 determinations of the same sample at the same day showed a Coefficient of Variation of 4.21%

Inter-Assay
The realization of 10 determinations of the same sample at different days showed a Coefficient of Variation of 4.46%

Analytical Specificity
A comparison with a reference method showed a correlation coefficient (r) of Q.989 obtained from a multi-atomic sample. The result equation of the linear regression is Y = 1,0041 X + 7,8468

CLINICAL SIGNIFICANCE

The plasma iron transport is carried out by a specific beta-globulin the Transferrin which may be determined by the iron quantity it can bind. This quantity is called the Total Iron Binding Capacity (TBC).

For the complete diagnostic of accumulation or deficiency of Iron in the patient, it's important to perform a complete haematological profile: iron determination, TBC and Transferrin Saturation Index (TSI).

The frequent reasons for the increasing TBC are: anaemia by iron deficiency, pregnancy, use of contraceptive pills, hypoxia events, etc. The TBC is decreased in chronic inflammatory diseases, malignant diseases, nephropathies, iron intoxication, talasemia major. The TSI is increased in iron intoxication, nephropathies, talasemia major, primary iron deficiency, haemolytic and sideroblastic anaemias. The TSI is decreased in iron deficiency, malignant diseases, chronic infections, pregnancy and menstrual periods.

OBSERVATIONS
1. The laboratory should establish chemical, microbiological and particulate requirements for water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be

B IOTÉCNICA INDCOM LTDA - CNPJ: 025340690001/20. Avenida Vashngton RBarço, 200 - Indústria | Nguêre de Luca CEP 37072-030 | Ver g nha- MG BRAS L T d f ax +55 35 3214 4646 www.bioteccai.ind.br

validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
2. The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

PRESENTATIONS 2
3

QUALITY ASSURANCE/CUSTOMER TECHNICAL SERVICE
Before being approved for use Biotécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiry date stated in the label of the external packaging, since it's stored and transported in the specified conditions.
The quality control data concerning this product or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@bioteccai.com.br

AUTOMATION
This product is compatible to the most types of biotechnical automatic analysers. The applications are available at: www.bioteccai.ind.br

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- GOODWIN J.; MURPHY, B.; GUILLETTE M Direct measurement of serum iron and binding capacity. *Clin Chem* v.12, p 47-57, 1966
- STOCKEY, L L. Ferroni - a new spectrophotometric reagent for iron. *Anal. Chem* v.42, p 779-781, 1970.
- ITANQMI CAP Comprehenhensiv e che mistry: serum iron survey. *Am J. Clin Pathol.* v.70, p516-522, 1978.
- YOUNG D.S Effects of drugs on clinical laboratory tests - vol. 2 5ed. Washington DC: AACC Press, 2000.
- WESTGARD J. O et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin Chem* v.27 p 493-501, 1981.

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