Cloretos

Chloride / Chloruros

Ref. 12.003.00

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Before being released for consumption, all Biotecnica reagents are tested by the Quality Control Department. The quality of the reagents is ensured by the expiration date mentioned on the presentation packaging since they are stored and transported under specified conditions. The data related to the Quality Control of this product (lot printed on the labels of the reagent bottles) or any doubt in the use of this kit, contact the Scientific Advisory Board of Biotécnica Ltda, by phone +55 35 3214 4646 or by email sac@biotecnicatda.com.br

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ENGLISH

INTENDED USE



Kit intended to determination of chloride in serum, plasma and urine. In vitro Diagnostic use only.

STORAGE AND HANDLING

- Store at room temperature.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

Chloride ions react with mercury thiocyanate to form mercuric chloride and free thiocyanate ions. They react with ferric ions to form ferric thiocyanate whose absorbance can be measured at 500 nm. The intensity of the color is proportional to the concentration of chlorides in the sample.

SAMPLES: TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum, plasma (heparin) and urine.

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material. Serum or plasma should be separated within 1 hour after collection. **Preservation:** Chloride in serum and plasma (heparin) is stable for 7 days at a temperature of 4 to 8 °C or 20 to 25 °C and 1 year at -20 °C. In urine is stable for 14 days at 2 to 8 °C or at room temperature and 60 days at - 20 °C.

Urine: Use sample collected within 24 hours, with no need to add preservatives. Mix well and measure the volume. Centrifuge an aliquot of the urine. Dilute an aliquot of urine supernatant centrifuged at 1: 2 ratio with purified water (where necessary, dilution of urine should be altered to give result within the operating range of the method). Use diluted urine for testing. Multiply the result obtained by 2 (or by the factor that has been applied).

PRODUCT DESCRIPTION

R 1	Mercury thiocyanate ≥ 1.0 mmol/L, Ferric nitrate ≥ 20 mmol/L, Nitric acid ≥ 0.1% v/v, sulfuric acid ≥ 0.1% v/v, mercury nitrate ≥ 0.001% v/v.	
STD	Chloride ions at a concentration equivalent to 100 mEq/L and preservative. Traceable to NIST certified reference material SRM 3182 and NIST SRM 999b.	

STABILITY IN USE

- The stability of the product (R1 and STD) in use is 24 months since the recommended storage conditions (room temperature) are followed.

PRODUCT TREATMENT AND HANDLING

A) REAGENT PREPARATION

R1 e STD: Ready-to-use reagents.

B) OPERATING RANGE

The operating range of the product is 66,87 mEq/L at 150,00 mEq/L.

For higher values, dilute the sample with purified water, carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For Calibration and Internal Quality Control Laboratory it is recommended to use the calibrator serum and the control sera below

Serum Calibrator - Autocal H	REF	13.002.00
Normal Control Serum - Quantinorm		13.003.00
Pathological Control Serum - Quantialt		13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the assay tubes:

	Blank	Standard	Sample
STD	-	5 µL	-
Sample	-	-	5 µL
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize firmly and keep the tubes for 5 minutes at room temperature.

3. Measure the absorbance of the Standard and the Sample against the Blank at 500 nm. The color is stable for 30 minutes.

B) CALCULATIONS

Chlorides (mEq/L) = $\frac{\text{Sample Absorbance}}{\text{Standard Absorbance}} \times \text{Standard Concentration (mEq/L)}$

Example:

Concentration of STD = 100 mEq/L

Sample Absorbance = 0,336

Absorbance of STD = 0,321

Chlorides (mEq / L) = $\frac{0,336}{0,321} \times 100 = 105 \text{ mEq/L}$

0,321

With Calibration Factor:

Calibration Factor = $\frac{\text{STD Concentration (mEq / L)}}{\text{Absorbance of STD}}$

Chlorides (mEq/L) = Sample Absorbance x Calibration Factor

Example: Calibration Factor = $\frac{100}{0,321} = 312$

Chlorides (mEq / L) = 0.336 x 312 = 105 mEq/L
24 Hour Urine Calculation

mEq/ 24 hours = * mEq/L x urinary volume of 24 h in liter

* Value corrected by dilution factor

C) INTERPRETATION

Chloride is the major extracellular anion. Together, sodium and chloride represent the bulk of the osmotically active constituents of plasma. The chloride ions from the feed are almost completely absorbed in the intestinal tract. It is significantly involved in: maintenance of water distribution, osmotic pressure and anion-cation balance in the extracellular fluid (ECF). Determination of plasma chloride concentration is useful in the differential diagnosis of acid-base disorders and is essential for calculating the anionic range. Chloride fluctuations in serum or plasma have little clinical consequence, but serve as signs of a basic disorder in fluid homeostasis and acid-base balance, and may assist in differentiating the cause of these disorders. Hyperchloremia occurs in dehydration, renal tubular acidosis (RTA), acute renal failure, metabolic acidosis associated with prolonged diarrhea and loss of sodium bicarbonate, diabetes insipidus, conditions of adrenocortical hyperfunction, overload of saline solutions and respiratory alkalosis as a result of renal excretion of HCO₃⁻. Hypochloremia is observed in metabolic acidosis caused by increased production or decreased excretion of organic acids (eg, diabetic ketoacidosis and renal insufficiency), persistent gastric secretion, prolonged vomiting, hypochloremic alkalosis, and total body chloride depletion with corresponding retention of HCO₃⁻. Other associated conditions include salt loss nephritis, Addison's disease, hyperaldosteronism, and bromide poisoning. The excess chloride is excreted in the urine and is also eliminated in sweat. Usually urinary excretion of chloride approaches dietary intake. Physiologically increased levels of post-menstrual diuresis and decreased pre-menstrual water retention of salt and water are observed.

INTERFERING AND LIMITATIONS

Hemolysis, jaundice and Lipemia: Hemoglobin > 500 mg/dL / Bilirubin > 40 mg/dL / Triglycerides > 1250 mg/dL interferes with the dosage.

Medications: Refer to recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Detection limit: 3,05 mEq/L / Limit of quantification: 66,87 mEq/L.

Analytical Specificity: The product specifically determines the chloride in the presence of other interfering substances in the sample, up to the concentrations reported above.

Correctness (Accuracy): The method was compared with a similar method by the determination of 40 samples in duplicate. The regression equation $y = 0.987x + 1.268$ and correlation coefficient $r=0.9972$ was obtained. Using this equation the total systematic error estimated was 0.51% for a level of 70.00 mEq/L and -0.24% for a level of 120.00 mEq/L.

Accuracy: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, obtaining:

Samples (mEq/L)	Repetitions	Within-Run Precision		Total Precision	
		SD (mEq/L)	%CV	SD (mEq/L)	%CV
86,1413	80	0,7691	0,9	0,8363	1,0
95,2131	80	0,8045	0,8	0,8523	0,9
135,9241	80	0,8150	0,6	0,8927	0,7

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with Good Clinical Laboratory Practices.
- Do not mix reagents from different lots.
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any sign in disagreement with the one specified in the product MSDS.
- Use glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.
- R1 is toxic because it contains Mercury thiocyanate.

REFERENCE RANGES

Serum and Plasma	98 - 107 mEq/L
Urine	110 - 250 mEq/24 hours

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

Conversion to International System Unit (SI):

Chlorides (mEq/L) X 1,0 = Chlorides (mmol/L)

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer able to read at 500 nm (420-520).
- Glass Pipettes and/or automatic.
- Clock or chronometer.
- Assay tubes.

WARNINGS AND PRECAUTIONS RELATED TO DISPOSAL OF THE PRODUCT

Disposal, Safety and First Aid information are described in the Individual Safety Data Sheet for this product, available at www.biotechnica.ind.br or by telephone +55 35 3214 4646.

Dispose of leftover reactions in accordance with Good Clinical Laboratory Practices (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use, BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotechnica.com.br

AUTOMATION

This procedure is automated in most analyzers. The protocols are available at www.biotechnica.com.br.

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