## Colinesterase

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#### ENGLISH

#### INTENDED USE

Kit for the quantitative determination of the enzymatic activity of Cholinesterase in serum and plasma. In vitro diagnostic use.

#### STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

#### WORKING PRINCIPLE

Cholinesterase catalyzes the hydrolysis of butyrylthiocholine forming butyrate and thiocholine. The Thiocholine reduces hexacyanoferrate III (yellow) to hexacyanoferrate II (colorless). The decrease in absorbance measured at 405 nm is directly proportional to the catalytic activity of Cholinesterase in the sample.

#### SAMPLES – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

**Sample Type:** Serum and plasma (EDTA and heparin).

**Collection, handling and preparation:** Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

**Preservation:** Cholinesterase in serum and plasma is stable for seven days if stored at 4 to 8 °C and one year if stored at -20 °C.

#### PRODUCT DESCRIPTION

<b>R 1</b>	Sodium pyrophosphate $\geq$ 20 mmol / L; Potassium ferricyanide $\geq$ 0,5 mmol / L and stabilizer.
------------	---

<b>R 2</b>	S-butryrylcholine $\geq$ 20 mmol / L and stabilizers.
------------	---

Traceable to specific absorptivity of chromogen hexacyanoferrate III using calibrated pipettes and handheld spectrophotometer that provide absolute values.

#### STABILITY IN USE

- After opening, the product (R1 and R2) in use is stable to the expiration date printed on the label, since recommended storage conditions (2 to 8 °C) are followed.
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

#### TREATMENT AND HANDLING OF THE PRODUCT

##### A) REAGENT PREPARATION

R1 and R2: Ready-to-use reagents.

##### B) OPERATING INTERVAL

The operating range of the product is 348.08 U / L a 15000 U / L.

For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

#### QUALITY CONTROL

The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For internal control of laboratory quality it is recommended to use the control sera below:

Normal Control Serum- Quantinorm	13.003.00
Serum Pathological Control - Quantialt	13.004.00

#### TEST PROCEDURE, CALCULATION AND INTERPRETATION

##### A) TEST PROCEDURE

1. Assay the samples and blank (purified water).
2. Pipette into test tube.

Purified water/sample	20 µL
R1	1000 µL

3. Homogenize and incubate at 37 °C for 3 minutes. Add:

R2	200 µL
----	--------

4. Homogenize and insert into the cuvette holder thermostated at 37 °C, zeroed with purified water at 405nm. Set the stopwatch.
5. After 90 seconds, note the initial absorbance A<sub>0</sub> and perform new readings after 30, 60 and 90 seconds (A<sub>1</sub>, A<sub>2</sub> and A<sub>3</sub> respectively).

##### B) CALCULATIONS

Using the absorbance readings, determine the ΔA / min of the sample and blank as follows:

$$\Delta A / \text{min} = \frac{(A_0 - A_1) + (A_1 - A_2) + (A_2 - A_3) \times 2}{3}$$

Cholinesterase activity in the sample is calculated with the following equation:

$$\text{Cholinesterase (U / L)} = [(\Delta A / \text{min of Sample} - \Delta A / \text{min of White})] \times 65800$$

#### C) INTERPRETATION

Cholinesterase is an enzyme found in the liver, pancreas, heart, serum and the white matter of the brain, but its physiological role is unknown. Measurement of cholinesterase in serum is useful as an indicator of possible insecticide intoxication, for the detection of patients with atypical form of the enzyme, or as a liver function test. Decreased levels of cholinesterase are found in cases of organophosphorus poisoning, cirrhosis, myocardial infarction, acute infections and atypical phenotypes of the enzyme.

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**Analytical Specificity:** The product specifically determines cholinesterase in the presence of other interfering substances in the sample, up to the concentrations reported above.

**Accuracy:** The method was compared with a similar method by the determination of 40 samples in duplicate. The regression equation Y = 1,000x-6,190 and correlation coefficient r = 0.9944 was obtained. Using this equation the estimated total systematic error of -0.100% for a level of 6000 U / L and -0.044% for a level of 14000 U / L.

**Accuracy:** It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, being obtained:

Samples (U/L)	Repetitions	Intra-run precision		Total Precision	
		SD (U/L)	%CV	SD (U/L)	%CV
1520,40	80	7,30	0,5	27,79	1,8
4101,98	80	71,458	1,7	76,446	1,9
7163,69	80	37,637	0,5	45,814	0,6

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

#### RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the

required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.

- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

#### REFERENCE RANGES

37 °C
4900 – 11900 U/L

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

Conversion to International System Unit (SI):  $\mu\text{kat} / \text{L}$   
Cholinesterase (U / L) x 0.0167 = Cholinesterase ( $\mu\text{kat} / \text{L}$ )

#### MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer for reading at 405 nm.
- Water bath, thermostated at 37 ° C.
- Glass and / or automatic pipettes.
- Clock or Stopwatch.
- Test tubes

#### ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: [sac@biotecnicaltda.com.br](mailto:sac@biotecnicaltda.com.br)

#### AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at [www.biotecnicaltda.ind.br](http://www.biotecnicaltda.ind.br)

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