PERFORMANCE CHARACTERISTICS

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Analytical Specificity: The product determines Lactate Dehydrogenase specifically in the presence of other interfering substances in the sample, up to the concentrations stated above.

Accuracy: The method was compared with a similar method for the determination of 40 samples in duplicate. The regression equation obtained was $y = 0,999x - 6,54$ and the correlation coefficient $r = 0,9958$. Using this equation the total systematic error estimates is -1,73% to a level of 400 U/L and -0,64% to a level of 1200 U/L.

Precision: Determined using samples in 03 levels of decision, with two runs in duplicate per day during 20 days, being obtained:

Sample (U/L)	Repetitions	Within-Run Precision		Total Precision	
		SD	%CV	SD	%CV
261,059	80	2,709	1,0	2,834	1,1
776,150	80	2,239	0,3	2,496	0,3
1320,326	80	8,64	0,7	9,191	0,7

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with Good Laboratory Practice
- Do not mix reagents from different lots.
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any significant disagreement with the one specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions, when they are not in use.
- Use glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The level of the water bath must be greater than the test tubes containing the reaction.
- The level of water in the thermostat bath must be over the level of mix reaction with the assay tubes.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

REFERENCE RANGES

ADULTS		37 °C
Men		214 - 450 U/L
Women		195 - 453 U/L

CHILDREN		37 °C
0 a 6 months		613 - 1020 U/L

These values are intended for orientation only. It's recommended that each lab establishes its own reference range.

Conversion for International System of Units (SI):

$$LDH(U/L) \times 0,017 = LDH(\mu\text{kat/L})$$

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer with a thermostated cuvette able to read at 340 nm.
- Thermostatic water bath at 37°C
- Glass pipettes and/or automatic.
- Clock or chronometer.
- Assay tubes

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.bioteccalind.com or calling for 55-35-3214-4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use Biotecnia reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiry date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@bioteccaltda.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.bioteccaltda.com

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