SANGUE OCULTO

FOB (Fecal Occult Blood)/ Sangre Oculta
Ref. 30.011.00

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ENGLISH

INTENDED USE

Kit for the qualitative determination of human occult blood in stool. Use *in vitro* diagnostic.

STORAGE AND HANDLING CONDITIONS

- Store at 2 to 30 °C.
- The validity of the kit is printed on the package label.
- Do not use the reaction strip if expire date is over.

WORKING PRINCIPLE

Chromatographic immunoassay on nitrocellulose membrane. The human hemoglobin present in the sample binds to an antibody conjugated with colloidal gold and migrates through the membrane by capillarity to the test region containing human anti-hemoglobin antibodies forming a colored band. The unbound conjugate continues to migrate to the control region where a colored band should always appear.

SAMPLES: TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Stool.

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. Samples to be analyzed should be treated as potentially infectious material

Preservation: The sample is stable for 3 days at a temperature of 2 to 8 °C.

PRODUCT DESCRIPTION

Device Single-use nitrocellulose strip-containing casing coated with colloidal gold-conjugated monoclonal antibodies, anti-human hemoglobin monoclonal antibodies.

Diluent Stabilizers.

STABILITY IN USE

- After opening the individual packaging the product is stable for 1h under the following humidity / temperature conditions: 35% / 20 °C, 48% / 24,6 °C and 85% / 37 °C.

PRODUCT TREATMENT AND HANDLING

A) REAGENT PREPARATION

Device and Diluent: Ready for use.

B) OPERATING INTERVAL

Up to 4000 µg/mL.

QUALITY CONTROL

The device has a control line that must be visible for the interpretation of the result.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

Before opening the packaging of the test device and the diluent, leave them at room temperature for 30 minutes (15 to 30 °C).

Dilution of sample:

1. Insert the stick in at least 3 different locations in the sample.
2. Return the stick to the collection tube, close tightly and shake vigorously. Wait 2 minutes.

Performance of the test:


1. Remove the device from its packaging and leave it on a horizontal surface.
2. Unscrew the cap from the diluent and add 2 drops (80 µL) of the diluted sample into the application scoop (S). Set the stopwatch.
3. Observe the result at exactly 5 minutes

B) ANALYSIS OF RESULTS

POSITIVE: two colored lines appear (red or pink), one line in the control region (C) and another in the test region (T).

NEGATIVE: A colored line appears only in the control region (C).

INVALID: no colored line appears in the control region (C).

C) INTERPRETATION

Investigation of bleeding in the gastrointestinal tract through stool occult blood research is an aid to clinical diagnosis and monitoring of bleeding disorders of the digestive tract. It is an indicated analysis for the preventive and anticipated screening of colorectal diseases in asymptomatic patients.

INTERFERES OR LIMITATIONS

Specificity: No cross-reactions were observed with sheep, chicken, ox, horse and pig hemoglobin at concentrations less than 500 µg / mL.

Medications: consult recommended literature reference (Young, 2000).

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High dose Prozone effect: No high dose prozone effect at concentrations up to 4000 µg / mL human hemoglobin is observed.

Analytical Specificity: The product specifically determines human hemoglobin in the presence of interferents in the sample, up to the concentrations reported above.

Accuracy: To estimate the percentages of agreement, the method was compared with a similar one by the determination of 100 samples, obtaining:

Percentage of Positive Agreement: 96%

Negative Agreement Percentage: 100%

Precision: It was determined using a sample with a concentration of 0,05 µg / mL of human hemoglobin in 20 replicates. The results obtained were in agreement with those expected.

RESIDUAL RISKS, CARE AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not use the product when it has a visual characteristic that is not in accordance with the product's MSDS.
- Don't leave the product out of the specified storage conditions.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Only use test and buffer devices whose packages are perfectly sealed.
- The test reading should not exceed 5 minutes.

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Clock or Stopwatch.

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at www.biotechnica.ind.br or by telephone (35) -3214-4646.
- Disposal of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY / SAC GUARANTEE – CUSTOMER SERVICE

Before being released for consumption, all Biotécnica reagents are tested by the Quality Control Department. The quality of the reagents is ensured by the expiry date mentioned on the presentation packaging, since they are stored and transported under specified conditions. The data related to the Quality Control of this product (lot printed on the labels of the reagent bottles) or any doubt on the use of this kit, contact the Scientific Advisory Board of Biotécnica Ltda, by phone +55 35 3214 4646 or by email sac@biotechnicaltda.com.br.

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- MANDEL, Jack S. et al. Sensitivity, Specificity, and Positive Predictivity of the Hemocult Test in Screening for Colorectal Cancers. *Gastroenterology*, v. 97, p.597-600, 1989.
- CATOMERIS, Peter et al. Effect of Temperature and Time on Fecal Hemoglobin Stability in 5 Fecal Immunochemical Test Methods and One Guaiac Method. *Arch Pathol Lab Med*, v. 142, p.75-82, jan. 2018.
- RABENECK, Linda et al. Fecal immunochemical tests compared with guaiac fecal occult blood tests for population-based colorectal cancer screening. *Can J Gastroenterol*, v. 26, n. 3, p.131-147, mar. 2012.
- UCHIDA, Kazuo et al. Immunochemical detection of human blood in feces. *Clínica Química Acta*, v. 189, p.267-274, 1990.
- YOUNG, D.S. *Effects of drugs on clinical laboratory tests - vol. 2*, 5 ed. Washington DC: AAC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.

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