IgG Plus

IgG Plus / IgG Plus
Ref. 20.019.00

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ENGLISH

INTENDED USE

Kit intended for the determination of Immunoglobulin G (IgG) in serum and plasma. Use *in vitro* diagnostic.

STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- The validity of the kit is printed on the package label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

An immunoturbidimetric method in which anti-human immunoglobulin G antibodies form insoluble complexes when mixed with samples containing IgG. The turbidity intensity is proportional to the IgG concentration of the sample.

SAMPLES – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum and plasma (heparin or EDTA).

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: The immunoglobulin G in serum and plasma is stable for 8 months at a temperature of 4 to 8 °C or - 20 °C.

PRODUCT DESCRIPTION

R 1	Goat anti-human immunoglobulin G antibodies; Tris buffer pH 8.2; preservative.	
------------	--	---

STABILITY IN USE

- After opening, the product (R1) in use is stable until the expiration date on the label, since recommended storage conditions (2 to 8 °C) are followed.

- Reagents should remain outside the specified temperature only for the time required to perform the tests.

TREATMENT AND HANDLING OF THE PRODUCT

A) REAGENT PREPARATION

R1: Reagent ready for use.

B) OPERATING INTERVAL

The product's operating range is 444 mg / dL to 2640 mg / dL.

For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. For Calibration and Internal Laboratory Quality Control it is recommended to use the calibrator and control below:

Multicalibrador Plus	REF	21.011.00
Medium Multiparameter Control		21.003.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

REACTIONS PARAMETERS	
 TYPE OF REACTION	Final Point /Growing
 Wave length	540 nm
 Temperature	37 °C
 Sample volume/CAL*	2 µL
 R1* Volume	285 µL
 Incubation time	2 minutes
 Calibration	Calibration curve**

*The sample volumes / CAL and R1 can be modified as long as the established ratio is maintained.

**CALIBRATION CURVE:

Dilute the PLUS Multicalibrator in 9g / L NaCl as follows:

Dilution	1	2	3	4	5	6
 Calibrador (µL)	-	10	25	50	75	100
 NaCl 9g/L (µL)	100	90	75	50	25	-
 Fator	0	0,1	0,25	0,5	0,75	1

Multiply the IgG concentration indicated in the Multicalibrator Plus instruction by the corresponding factor to obtain the IgG concentration of each dilution. The IgG concentration of the calibrator is traceable to Standard Reference Materials ERM DA470k / IFCC.

B) CALCULATIONS

The equipment automatically calculates the analyte concentration in each sample.

C) INTERPRETATION

Quantification of immunoglobulins in serum is important in the diagnosis of primary or secondary immunodeficiency, in the monitoring of immunoglobulin therapy and in the monitoring of multiple myeloma. IgG accounts for about 75% of the total circulating immunoglobulins produced by plasma cells. The decrease in IgG levels occurs mainly due to congenital and acquired immunodeficiencies. Increases in IgG concentration are a consequence of polyclonal hyperglobulinemia, cirrhosis and autoimmune diseases. Increases in monoclonal IgG concentration are observed in multiple myeloma, lymphocytic leukemia and Waldenström macroglobulinemia.

INTERFERING AND LIMITATIONS

Hemolysis, Jaundice and Lipemia: Hemoglobin> 400 mg / dL, Bilirubin> 40 mg / dL and Triglycerides> 500 mg / dL, interfere with the assay.

Medications: see recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Detection limit: 1.48 mg / dL / Limit of quantification: 444 mg / dL.

Prozone High Dose Effect: Prozone effect does not occur up to 7000 mg / dL IgG. **Analytical Specificity:** The product specifically determines Immunoglobulin G in the presence of other interfering substances in the sample up to the concentrations reported above.

Accuracy: The method was compared with a similar method by the determination of 31 samples in duplicate. The regression equation $Y = 1,000X + 5,924$ and correlation coefficient $r = 0,9906$ was obtained. Using this equation the estimated total systematic error of 0.66% for a level of 900 mg / dL and 0.46% for a level of 1300 mg / dL.

Accuracy: It was determined using samples at 03 decimal levels, two runs in duplicate per day for 20 days, being obtained.

Samples (mg/dL)	Repetitions	In-run Precision		Total Precision	
		SD (mg/dL)	%CV	SD (mg/dL)	%CV
670	80	6,25	0,9%	8,18	1,2%
1245	80	8,83	0,7%	21,47	1,7%
2789	80	38,43	1,4%	41,34	1,5%

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the reagents out of the specified storage conditions when they are not in use.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

REFERENCE RANGES

Newborn (4 days)	700 to 1480 mg/dL
Adults (20-60 years)	700 to 1600 mg/dL
Adults (> 60 years)	600 to 1560 mg/dL

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer for reading at 540 nm.
- Glass and / or automatic pipettes.
- Test tubes
- Biotécnica Plus Multicalibrator – Ref. BT 21.011.00

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at www.biotecnica.ind.br or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotecnicaltda.com.br

AUTOMATION

This procedure is automated in most analyzers. The protocols are available at www.biotecnica.ind.br

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- BURTIS, C.A.; ASHWOOD, E.R.; BRUNS, D.E. Tietz: Fundamentos de Química Clínica. 6. ed. Rio de Janeiro: Saunders Elsevier, 2008. 959 p.
- DATI, F. *et al.* Consensus of a Group of Professional Societies and Diagnostic Companies on Guidelines for Interim Reference Ranges for 14 Proteins in Serum Based on the Standardization against the IFCC/BCR/CAP Reference Material (CRM 470). *Eur J Chem Clin Biochem.* v.34, p.517-520, 1996.
- MARIA, R.M; SILVA, R.O.P. Gamopatias monoclonais – critérios diagnósticos e diagnósticos diferenciais. *Rev. bras. hematol. Hemoter.* v.29, n.1, p.17-22, 2007.
- SILVA, M.T.N. *et al.* Imunoglobulinas séricas em crianças com exposição perinatal ao vírus da imunodeficiência humana. *Jornal de Pediatria.* v.77, n.3, p.209-218, 2001.
- SKOUG, J.W, PARDUE, H.L. Kinect Turbidimetric Method for the immunochemical Quantification of Immunoglobulins, Including Samples with Excess Antigen. *Clin. Chem.* v. 34, n. 2, p. 309-315, 1998.
- WESTGARD, J. O. *et al.* A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: **AACC Press**, 2000.

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