Potássio Enzimático

Enzymatic Potassium / Potasio Enzimático
Ref. 12.012.00

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ENGLISH

INTENDED USE

Kit for the determination of potassium in serum. In vitro diagnostic use.

STORAGE AND HANDLING CONDITIONS

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

Potassium is determined enzymatically through the activity of potassium-dependent pyruvate kinase, which reacts with the substrate phosphoenolpyruvate to produce pyruvate. This reacts with NADH in the presence of LDH to form Lactate and NAD. The disappearance rate of NADH at 340 nm is proportional to the potassium concentration in the sample.

SAMPLES – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum.

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: Potassium in the serum, is stable for 6 weeks stored at a temperature of 4 to 8 °C and 1 year at a temperature of -20 °C.

PRODUCT DESCRIPTION

R 1	Tris buffer 250 mmol / L pH 8.2; Screener 12 mmol / L; PEP ≥ 3.3 mmol / L; ADP ≥ 3.15 mmol / L; α-oxoglutarate ≥ 1.2 mmol / L; NADH ≥ 0.35 mmol / L; GLDH ≥ 11 U / mL; PK ≥ 1.2 U / mL.	
R 2	LDH ≥ 65 U/mL.	
STD 1	KCl solution with potassium concentration indicated on bottle label	
STD 2	KCl solution with potassium concentration indicated on bottle label	

Traceability of the STD1 and STD2 standards was obtained with the GBW09152 secondary standard of concentration defined by ID-MS.

STABILITY IN USE

- After opening, the product (R1, R2, STD1, STD2) in use is stable to the expiration date printed on the label, since recommended storage conditions are followed (2 to 8 °C).
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

TREATMENT AND HANDLING OF THE PRODUCT

A) REAGENT PREPARATION

R1, R2, STD1 and STD2: Ready-to-use reagents.

B) OPERATING INTERVAL

The operating range of the product is from 0.37 mmol / L to 10 mmol / L.

For higher values, dilute the sample with purified water, carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For Internal Quality Control Laboratory it is recommended to use the controls below:

Normal Control Serum- Quantinorm **REF** 13.003.00
Serum Pathological Control - Quantial 13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1 - Pipette into test tube.

	STD1	STD2	Sample
R1	720 µL	720 µL	720 µL
Sample	20 µL	20 µL	20 µL

2. Homogenize and incubate at 37 °C for 5 minutes.

3. Add:

R2	240 µL	240 µL	240 µL
----	--------	--------	--------

4. Homogenize and transfer immediately to a thermostated cuvette. Set the stopwatch. Note the absorbance at 120 seconds (A1) and 240 seconds (A2) at 340 nm by zeroing with purified water.

B) CALCULATIONS

Potassium (mmol/L) = (Δ Sample Absorbance x Calibration Factor) – Intersection

$$\text{Calibration Factor} = \frac{[\text{STD2}] - [\text{STD1}]}{\Delta\text{Abs STD2} - \Delta\text{Abs STD1}}$$

Intersection = Calibration Factor x Δ Abs STD 1 - [STD1]

Where:

Δ Sample absorbance = (A₁ – A₂) from sample

[STD2] = STD2 concentration

[STD1] = STD1 concentration

Δ Abs STD2 = (A₁ – A₂) from STD2

Δ Abs STD1 = (A₁ – A₂) from STD1

Exemplo:

	Absorbance		ΔAbs
	Abs 1	Abs 2	A1-A2
STD1	3,0 mmol/L	1,670	1,423
STD2	7,0 mmol/L	1,549	1,168
Sample		2,173	1,887

$$\text{Calibration Factor} = \frac{7,0 - 3,0}{0,381 - 0,247} = \frac{4,0}{0,134} = 29,85$$

Intersection = (29,85 x 0,247) – 3,0 = 4,37

Sample concentration (mmol/L) = (0,286 x 29,85) – 4,37 = 4,17

C) INTERPRETATION

The function of electrolytes in the human body is multiple. These include maintenance of osmotic pressure and water distribution in the body's various fluid compartments, maintenance of adequate pH, regulation of heart and other muscle function, participation in oxidation-reduction reactions, and catalysis as enzymatic cofactors. Abnormal levels of electrolytes may be the cause or the consequence of various disorders. Potassium is the major cation of intracellular fluid. Its intracellular concentration is maintained by the Na-K-ATPase pump, which continuously carries potassium into the cell. This pump is responsible for the adjustments of the ionic gradients on which the nerve impulses and contractility of the muscle depend. Elevated levels of potassium (hyperkalemia) are present in different clinical conditions such as: metabolic acidosis, dehydration, obstructive nephropathy, use of potassium-sparing diuretics, etc. Low levels of potassium (hypokalemia) occur in alkalosis / metabolic acidosis, extra-renal loss (diarrhea, excessive sweating), poor diet in potassium, among others.

INTERFERING AND LIMITATIONS

Hemolysis, Jaundice and Lipemia: Hemoglobin > 200 mg / dL / Bilirubin > 40 mg / dL / Triglycerides > 1250 mg / dL.

Medications: consult recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Limit of detection: 0.22 mmol / L / Limit of quantification: 0.37 mmol / L.

Analytical Specificity: The product specifically determines Potassium in the presence of other interfering substances in the sample up to the concentrations reported above.

Accuracy: The method was compared with a similar method by the determination of 40 samples in duplicate. The regression equation $y = 0.997x - 0.001$ and correlation coefficient $r = 0.9905$ was obtained. Using this equation the estimated total systematic error of -0.33% for a level of 3.80 mmol / L and -0.31% for a level of 8.00 mmol / L.

Accuracy: it was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, being obtained.

Samples (mmol/L)	Repetitions	In-run Precision		Total Precision	
		SD (mmol/L)	%CV	SD (mmol/L)	%CV
3,67	80	0,05	1,5	0,06	1,6
5,71	80	0,08	1,4	0,09	1,5
6,83	80	0,10	1,4	0,10	1,5

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it exhibits a visual characteristic in disagreement with that specified in the product's MSDS.
- Avoid leaving reagents out of the specified storage conditions when they are not in use.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator, and reagent.
- The water level of the water bath should be higher than that of the test tubes containing the reactions.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the requirements by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

REFERENCE RANGES

Serum: 3.5 to 5.1 mmol / L

Conversion: SI Unit (mmol / L) x 1 = Conventional Unit (mEq / L).

These values are for guidance only and it is recommended that each laboratory establishes its own reference range).

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer for reading at 340 nm
- Glass and / or automatic pipettes.
- Clock or Stopwatch.
- Test tubes

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at www.biotechnica.ind.br or by telephone (35) 3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotechnicaltda.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.biotechnicaltda.ind.br

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